

(0086571-86822701



Add.: # 701, Building C, Wanda Plaza, No. 666,
Hangxing Road, Gongshu District,
Hangzhou, 310015, Zhejiang Province, China

Web: www.weheremedical.com


Wehere

Hangzhou Wehere Medical Technology Co., Ltd.

Chest Drainage Tube



Drainage Tube - Chest Type



Basically designed for quick non operative chest as well as chest drainage.

Very soft tangle resistance tubing allows equal flow rate.

Cross side eyes to avoid aspiration of tissue.

Wide smooth drainage eyes ensure efficient drainage.

Proximal end is equipped with tapered connector for convenient connection for draining bottle.

X-ray opaque line given throughout.

Tube Without Trocar

Open distal end with 2 lateral eyes
allow non-traumatic insertion.

Size: FR 08- 32

PVC material (REF 0304A) or silicone
material (REF 0304B)are available.

Tube With Trocar

open distal end with 2 lateral eyes
ensure non-traumatic insertion.

Size: FR 08- 36

PVC Material (REF 0304C)