



INTRODUCTION

Hangzhou Wehere Medical Technology Co., Ltd. is a company on the medical health products manufacturing and trading in China. We are professional medical supplies nearly twenty years. On the basis of production management MDD and ISO13485 in our owner factories, guided by brand strategy, under the direct participation of American companies, we are committed to becoming an important member of the global medical and health-care industry.

We are here to support your medical & healthcare needs. Your any question are welcome!


Vision

To become the most respected, competitive and influencial medical device enterprise in China and give more people access to better health services.

**WEHERE
MEDICAL**
we here, we care



急救创伤绷带 (EMERGENCY BANDAGE)



HF FAB series

- Materials: Bandage + Cotton pad
- Color: Various
- Cotton pad size: 10cmx18cm; 15cmx18cm etc
- Width: 10cm, 15cm, 20cm
- Length (stretched): 3.6m, 4m, 4.5 etc
- Function: It is mainly used in hospital operating room/feld/road traffic to rescue the wounded.

Green Color & Two Pad Type



Skin Crepe & Hook



Skin Crepe & Bowl



Cohesive



Skin Color & Hook Strip



Horn & Bowl



Black High Elastic



医用绷带 MEDICAL BANDAGE

皱纹绷带 CREPE BANDAGE

- Materials: 80% cotton; 20% spandex
- Weight gsm: 60g/m², 65g/m², 70g/m², 75g/m², 80g/m² etc.
- Specification: Width 5cm, 7.5cm, 10cm, 15cm, 20cm
Length (stretched) 3.5m, 4m, 4.5m, 5m etc.
- Clip: With or without clips, elastic band clips, metal clips
- Basic packing: Individually packed in cellophane.
- Function: Comfortable to wear, strong air permeability, protection/fixation



HFA-1 氨纶本白皱纹绷带
Natural color elastic crepe bandage



HFA-2 氨纶肤色皱纹绷带
Skin color elastic crepe bandage



HFA-3 氨纶漂白皱纹绷带
Bleached elastic crepe bandage



HFA-4 氨纶本白红蓝线皱纹绷带
Red (Blue) line elastic crepe bandage



HFA-5 氨纶特白蓝线皱纹绷带
Bleached blue line elastic crepe bandage



HFB-1 氨纶本白平纹绷带
Natural plain elastic bandage



HFB-2 氨纶肤色平纹绷带
Skin color plain elastic bandage



HFB-3 氨纶漂白平纹绷带
Bleached plain elastic bandage

平纹绷带 PLAIN BANDAGE

- Materials: 55% cotton; 45% spandex
- Weight gsm: 70g/m², 75g/m², 80g/m², 85g/m² etc.
- Specification: Width 5cm, 7.5cm, 10cm, 15cm, 20cm
Length (stretched) 3.5m, 4m, 4.5m, 5m etc.
- Clip: With or without clips, elastic band clips, metal clips
- Basic packing: Individually packed in cellophane.
- Function: Comfortable to wear, strong air permeability, protection/fixation

医用绷带 MEDICAL BANDAGE

PBT绷带 PBT (CONFORMING) BANDAGE

- Materials: 100% polyester / 55% cotton; 45% polyester
- Weight gsm: 28g/m², 30g/m², 55g/m² etc.
- Specification: Width 5cm, 7.5cm, 10cm, 15cm, 20cm
Length (stretched) 3.5m, 4m, 4.5m, 5m etc.
- Clip: With or without clips, elastic band clips, metal clips
- Basic packing: Individually packed in cellophane.
- Function: Comfortable to wear, strong air permeability, protection/fixation



HFC-1 漂白PBT绷带
Bleached PBT bandage



HFC-2 皱纹PBT绷带
Crepe PBT bandage



HFC-3 厚型PBT绷带
Thick PBT bandage



HFC-4 PBT急救绷带
PBT First aid bandage



HFC-5 自粘皱纹PBT绷带
Cohesive crepe PBT bandage



高弹绷带 High Elastic Bandage



HF A1 肤色橡筋高弹绷带
Skin color high elastic bandage
● Materials:
70% polyester; 30% rubber/ 90% polyester;
10% spandex



HF B1 肤色高弹绷带
Skin color high elastic bandage
● Materials:
85% polyester;
15% rubber



HF F1 肤色高弹绷带
Skin color high elastic bandage
● Materials:
57% cotton;
26.5% polyester; 16.5% rubber



HF K1 开司棉高弹绷带
Purpled high elastic bandage
● Materials:
55% KS cotton;
20% polyester; 25% rubber



HF N1 花边高弹绷带
Laced high elastic bandage
● Materials:
50% cotton;
10% spandex; 40% polyester



HF B2 高弹绷带
High elastic bandage
● Materials:
75% polyester; 25% rubber / 90% polyester;
10% spandex



HF B3 蓝线高弹绷带
Blue line high elastic bandage
● Materials:
50% cotton;
20% polyester; 30% rubber




HF H1 高弹绷带
High elastic bandage
● Materials:
50% cotton;
20% polyester; 30% rubber

自粘绷带 COHESIVE (ELASTIC) BANDAGE



HF ECB-1 自粘绷带 (无纺布/全棉) Cohesive bandage (Non-woven/Cotton)

- Materials: Materials: Elastic fiber with non-woven/cotton
- Color: Various
- Specification: Width 2.5cm, 5cm, 7.5cm, 10cm, 15cm
Length (stretched) 4m, 4.5m, 5m etc.
- Basic packing: Individually packed in cellophane
- Function: Comfortable to wear, strong air permeability, protection/fixation, high self-adhesive



HF KT-1 运动胶带 Kinesiology Tape

- Materials: Cotton/Synthetic
- Color: Various
- Specification: Width 2.5cm, 3.8cm, 5cm, 7.5cm, 10cm; Length 5m
- Basic packing: Individually packed in cellophane
- Function: Comfortable to wear, strong air permeability, protection/fixation, high self-adhesive



HF EAB-1 EAB粘性绷带
Elastic Adhesive bandage (EAB)

- Materials: Cotton
- Kind of glue: Hot melt glue with Zinc Oxide
- Specification: Width 2.5cm, 5cm, 7.5cm, 10cm, 15cm
Length (stretched) 4m, 4.5m, 5m etc.
- Color: White or Yellow
- Function: Comfortable to wear, strong air permeability, protection/fixation, high self-adhesive

管状绷带 TUBULAR GAUZE BANDAGE

- Materials: 100% cotton / 100% polyester
- Weight gsm: 200g/m²
- Color: White or Yellow
- Function: Comfortable to wear, strong air permeability, protection/fixation, high self-adhesive



HFT-1 Tubular Gauze Bandage / Stockinette (100% cotton)

Product	Item No.	Width	Recommended use (body part)
Tubular gauze bandage 100% cotton	HFT-101	5cm	Small limbs
	HFT-102	7.5cm	Medium / Large limbs
	HFT-103	10cm	Extra large limbs
	HFT-104	15cm	Adult trunks

Noted: Sizes and color boxes can be customized

HFT-2 Tubular Gauze Bandage (100% polyester)

Product	Item No.	Width	Recommended use (body part)
Tubular gauze bandage 100% cotton	HFT-201	5cm	Small limbs
	HFT-202	7.5cm	Medium / Large limbs
	HFT-203	10cm	Extra large limbs
	HFT-204	15cm	Adult trunks
Noted: Sizes and color boxes can be customized			

HFT-3 Elastic Tubular Gauze Bandage (Latex)

Function:

- Provide comfortable and effective organizational
- Reduces patient pain and promotes healing
- Easy to use, strong flexibility
- Wide range of use, suitable for all parts of the human body
- Not brushed, not broken



Product	Item No.	Width	Recommended use (body part)
Elastic tubular gauze bandage	HFT-301	4.5cm	Small limbs
	HFT-302	6.25cm	Hands / Feet / Legs
	HFT-303	6.75cm	Arms / Legs
	HFT-304	7.5cm	Knees / Thighs
	HFT-305	8.75cm	Knees / Thighs
	HFT-306	10cm	Thighs
	HFT-307	12cm	S Trunks
	HFT-308	17.5cm	M Trunks
	HFT-309	21.5cm	L Trunks
	HFT-310	32.5cm	S Limbs
Noted: Sizes and color boxes can be customized			

网状绷带 TUBULAR NET BANDAGE

Item No.	Size	Image	Pcs/ctn	Box size	Used position
HFT-4 #0.5	0.8CM×25M		180	18*7*35cm	Finger
HFT-4 #1	1.7M×25M		120	18*13*2cm	Finger Legs Wrist
HFT-4 #2	2CM×25M		120	18*13*2cm	Hand Wrist
HFT-4 #3	2.3CM×25M		84	18*13*3cm	Arm Elbow
HFT-4 #4	2.7CM×25M		84	18*13*3cm	Foot Neck
HFT-4 #5	3CM×25M		84	18*13*3cm	Knee Legs
HFT-4 #6 #7	3.5CM×25M 4.2CM×25M		56	18*13*4.5cm	Head Legs Shoulder Chest
HFT-4 #8 #9 #10	5.8CM×25M 7.2CM×25M 8.5CM×25M		32	18*13*8cm	Buns Backside Buns
HFT-4 #11 #12	9CM×25M 10CM×25M		23	18*13*10cm	Large backside Chest

石膏绷带 PLASTER OF PARIS BANDAGE (POP)

HF P-1 石膏绷带 Plaster of Paris Bandage (POP)

- Materials: Cotton and plaster of paris
- Color: white
- Specifcation: Width 5cm, 7.5cm, 10cm, 15cm, 20cm
Length (stretched) 2.7m, 3m, 4m, 4.5m etc.
- Solidifying time: 2 to 5 minutes
- Basic packing: Individually packed in waterproof bag
- Function: Fixation of all kinds of fracture injuries, orthopedic and surgical procedures



石膏衬垫 ORTHOPEDIC BANDAGE (CAST PADDING)



Polyester



Cotton



Viscose

HF P-2 石膏衬垫 Orthopedic Bandage (Cast Padding)

- Materials: 100% polyester / 100% cotton / 50% viscose; 50% polyester
- Weight gsm: 75g, 80g, 90g etc.
- Specifcation: Width 5cm, 7.5cm, 10cm, 15cm, 20cm
Length (stretched) 2.7m, 3m, 4m, 4.5m etc.
- Basic packing: Individually packed in waterproof/paper bag
- Advantage: Preventing the plaster from coming into direct contact with the injured skin, protect skin from burns

纱布 COTTON GAUZE



HF D1 织边纱布绷带 Gauze Bandage with Woven Sides

- Material: 100% cotton
- Yarn: 40's, 32's, 21's
- Mesh: 19x15, 30x20, 30x30, 30x28, 40x20, 28x24
- Size: Length(stretched): 4m, 4.5m, 5m, 6m, 9m, 10m
- Width: 4cm, 5cm, 6cm, 7.5cm, 8cm, 10cm, 15cm
- Basic packing: Individually packed in cellophane
- Note: Personalized specification as possible at customer's Request.



HF Z-3 凡士林纱布片 Vaseline Gauze

- 10cmx10cm,
- 1pc/Sterile Package; 10packs/box
- Gamma Sterilization
- Vaseline Bandage
- 10cmx10m; 10cmx8m
- 1roll/Sterile Package; Gamma Sterilization



HF CZ01 无菌纱布片 Sterile/Non Sterile Gauze Swab

- Material: 100% cotton
- Yarn: 40's, 32's, 21's
- Mesh: 19x15, 20x12, 30x20, 24x20, 30x28, 12x8 etc
- Size: 2"x2"x8ply, 3"x3"x8ply, 4"x4"x8ply;
2"x2"x12ply, 3"x3"x12ply, 4"x4"x12ply;
2"x2"x16ply, 3"x3"x16ply, 4"x4"x16ply etc.
- Client request size is also workable.
- Package: 1pc/bag, 2pcs/bag, 100pcs/paperpack etc

纱布 COTTON GAUZE



HF CZ02 切边纱布绷带 Absorbent Gauze Bandage (WOW)

- Material: 100% cotton
- Yarn: 40's 32's 21's
- Mesh: 19x15,20x12,30x20,24x20,30x28,12x8 etc
- Size: Length:3.6m,4m,4.5m,5m,9.1m, 10m,13.5m etc
- Width: 2.5cm,4cm,5cm,7.5cm,8cm,10cm,.15cm,20cm etc
- Package: 1pc/paper roll.12rolls/pack.



HF CZ03 纱布卷 Gauze Roll

- Material:100% cotton
- Yarn: 40's 32's 21's
- X-Ray: X-ray detectable or not
- Mesh: 19x15,20x12,30x20,24x20,30x28,12x8 etc
- Size: 36"x5m,35"x6yards ,36"x50m,36"x100yards, 36"x100m etc
- Client request size is also workable.
- Package: 1pc/paper pack.



HF CZ04 手术巾 Laparotomy Sponge

- Material: 100% cotton
- Yarn: 40's 32's 21's
- X-Ray: X-ray detectable or not
- Mesh: 19x15,20x12,30x20,24x20,30x28,12x8 etc
- Size: 450x370mmx4ply/6ply, 350x350mmx4ply/6ply 200x170mmx4ply/6ply,350x370mmx4ply/6ply
- Client request size is also workable.
- Package: 1pc/paperpack.

医用胶带 MEDICAL TAPE / WOUND DRESSING

Non-Woven Tape(Skin/white)

Item No.	Size	Package	MOQ
HF Z-1	1.25cm × 5m	24rolls/box,30boxes/ctn	1500boxes
HF Z-2	2.5cm × 5m	12rolls/box,30boxes/ctn	
HF Z-3	5cm × 5m	6rolls/box,30boxes/ctn	
HF Z-4	7.5cm × 5m	6rolls/box,30boxes/ctn	
HF Z-5	10cm × 5m	6rolls/box,30boxes/ctn	

Transparent PE Tape

Item No.	Size	Package	MOQ
HF Z-6	1.25cm × 5m	24rolls/box,30boxes/ctn	1500boxes
HF Z-7	2.5cm × 5m	12rolls/box,30boxes/ctn	
HF Z-8	5cm × 5m	6rolls/box,30boxes/ctn	
HF Z-9	7.5cm × 5m	6rolls/box,30boxes/ctn	
HF Z-10	10cm × 5m	6rolls/box,30boxes/ctn	



医用胶带 MEDICAL TAPE / WOUND DRESSING

Silk Tape

Item No.	Size	Package	MOQ
HF Z-11	1.25cm × 5m	24rolls/box,30boxes/ctn	1500boxes
HF Z-12	2.5cm × 5m	12rolls/box,30boxes/ctn	
HF Z-13	5cm × 5m	6rolls/box,30boxes/ctn	
HF Z-14	7.5cm × 5m	6rolls/box,30boxes/ctn	
HF Z-15	10cm × 5m	6rolls/box,30boxes/ctn	

Zinc Oxide Adhesive Tape(Skin/White)

Item No.	Size	Package	MOQ
HF Z-16	1.25cm × 5m	24rolls/box,30boxes/ctn	1500boxes
HF Z-17	2.5cm × 5m	12rolls/box,30boxes/ctn	
HF Z-18	5cm × 5m	6rolls/box,30boxes/ctn	
HF Z-19	7.5cm × 5m	6rolls/box,30boxes/ctn	
HF Z-20	10cm × 5m	6rolls/box,30boxes/ctn	

PU/Non-Woven Fixing Tape

Item No.	Size	Package	MOQ
HF Z-21	5cm × 10m	1roll/box,160boxes/ctn	1500boxes
HF Z-22	10cm × 10m	1roll/box,80boxes/ctn	
HF Z-23	15cm × 10m	1roll/box,60boxes/ctn	
HF Z-24	20cm × 10m	1roll/box,40boxes/ctn	
HF Z-25	25cm × 10m	1roll/box,30boxes/ctn	
HF Z-26	30cm × 10m	1roll/box,25boxes/ctn	

医用胶带 MEDICAL TAPE / WOUND DRESSING

Sterile Wound Dressing(Non-Woven/PU)



Item No.	Size	Package	Item No.	Size	Package	Item No.	Size	Package
HF Z-27	5cm x 5cm	50pcs x 24boxes/carton	HF Z-37	6cm x 10cm	50pcs x 24boxes/carton	HF Z-47	9cm x 30cm	50pcs x 12boxes/carton
HF Z-28	5cm x 7cm	50pcs x 24boxes/carton	HF Z-38	7cm x 8cm	50pcs x 24boxes/carton	HF Z-48	9cm x 35cm	50pcs x 12boxes/carton
HF Z-29	5cm x 8cm	50pcs x 24boxes/carton	HF Z-39	7cm x 10cm	50pcs x 24boxes/carton	HF Z-49	10cm x 10cm	50pcs x 24boxes/carton
HF Z-30	5cm x 8.3cm	50pcs x 24boxes/carton	HF Z-40	8cm x 10cm	50pcs x 24boxes/carton	HF Z-50	10cm x 12cm	50pcs x 24boxes/carton
HF Z-31	5cm x 9cm	50pcs x 24boxes/carton	HF Z-41	8cm x 15cm	50pcs x 24boxes/carton	HF Z-51	10cm x 15cm	50pcs x 24boxes/carton
HF Z-32	5cm x 10cm	50pcs x 24boxes/carton	HF Z-42	9cm x 10cm	50pcs x 24boxes/carton	HF Z-52	10cm x 20cm	50pcs x 24boxes/carton
HF Z-33	5cm x 15cm	50pcs x 24boxes/carton	HF Z-43	9cm x 15cm	50pcs x 24boxes/carton	HF Z-53	10cm x 23cm	50pcs x 18boxes/carton
HF Z-34	6cm x 7cm	50pcs x 24boxes/carton	HF Z-44	9cm x 20cm	50pcs x 24boxes/carton	HF Z-54	10cm x 25cm	50pcs x 18boxes/carton
HF Z-35	6cm x 8cm	50pcs x 24boxes/carton	HF Z-45	9cm x 25cm	50pcs x 18boxes/carton	HF Z-55	10cm x 30cm	50pcs x 12boxes/carton
HF Z-36	6cm x 9cm	50pcs x 24boxes/carton	HF Z-46	9cm x 25cm	50pcs x 18boxes/carton	HF Z-56	15cm x 15cm	50pcs x 24boxes/carton

个人护理 / 防护 PERSONAL CARE / SAFETY

HFPC-1 颈托 Soft Collar

- Shaped to comfort and to give the correct support to relieve the pain suffered from whiplash, stiff neck and musclartension.
- Made from high density closed cell foam(4D&100D) for better support,have a woven covering.
- Available in 3size: small,medium and large
S: 44x9cm M:48x10cm L: 53x11cm



HFPC-2 吊带 Collar & Cuff

- Material: Cotton,Foam
- Size: 5cmx6m;7.5cmx6m
- Package: 2rolls/bag/box,16rolls/carton



HFPC-3 头颈固定器 Head&Neck Holder

- Material: Cotton, Foam
- Color: White, Skin
- Size: Normal



个人护理 / 防护 PERSONAL CARE / SAFETY

HFPC-4 皮肤牵引器 Skin Traction Kit

- Material: bandage and foam,belt and tape
- Size: child,Adult
- Type: Adhesive,Non adhesive
- Package: 1pc/bag/box,100boxes/carton



HFPC-5 鼻绷带 Nose Bandage

- Material: Nylon+Cotton
- Size: 2x45cm (Cotton size 8cmX1.8cm)
- Weight: 3.2±0.1g
- Elasticity: 180±20%
- Package: 20 or 40pcs/box



HF TB01 三角绷带 Triangular Bandage

- Material: 100% cotton,100% Non woven
- Size: 90cmx90cmx127cm/100cmx100cmx142cm
- Specification: Cotton 40S,32S,50x48,40x40 etc
- Nonwoven: 30g/m²,40g/m² etc



个人护理 / 防护 PERSONAL CARE / SAFETY

HF E01 冰绷带 Ice Bandage

- Pain relief in seconds
- Effective cold therapy for many conditions such as bruises, inflammations, sprains strains and sport injuries.
- No refrigeration necessary
- Hypoallergenic, good for hours, easy to use



HF N01 一次性医用短裤 Disposable Mesh Panty

- Material: Polyester and Spandex
- Use: Great for post-surgical & postpartum recovery, without adding pressure



HF AP01 一次性酒精片 Alcohol Pad

- Size: 30x60mm
- Material: Non-woven
- Color: White
- Packing: 50pcs/box, 100pcs/box
- Use: For disinfection purposes, such as skin disinfection, disinfection of articles, medical instruments, etc



个人护理 / 防护 PERSONAL CARE / SAFETY

HF FC-1 高分子绷带 Fiberglass Cast

- Material: Fiberglass, Polyester, PU foam
- Size: 2"/3"/4"/5"/6"x3.6m
- Package: 1pc/bag, 90rolls/120rolls/ctn
- Function: Light, breathable, high strength



HFS-1 口罩 Face Mask

- Material: Melt-blown non-woven PP fabric
- Color: Various
- Size: 17.5x9.5cm (Adult); 14.5x9.5cm (Kid) etc.
- BEF: 95%, 98%, 99%
- Packing: 50pcs/box, 40boxes/ctn



HFS-2 防护面罩 Medical Face Shield

- Material: PET
- PET Size: 32x22x0.25cm
- Packing: 1pc/bag, 240pcs/ctn



HFS-3 防护服 Isolation Gown (Coverall)

- Material: PP/SMS/SF/MP
- Size: S, M, L, XL, XXL, XXXL
- Weight: 35gsm, 40gsm, 50gsm, 60gsm etc.
- Packing: 1pc/pouch, 50pcs/ctn (sterile)
5pcs/pouch, 100pcs/ctn (non-sterile)

急救系列 (FIRST AID)

单兵急救包 IFAK/FIRST AID KIT

- | | |
|----------------------|------------------------------------|
| 1.Triangular bandage | 11.Alcohol pad |
| 2.Medical tape | 12.Elastic bandage |
| 3.Adhesive bandage | 13.Ax safety hammer |
| 4.Splint | 4.Outdoor folding shovel |
| 5.Whistle | 5.Multi-purpose knife |
| 6.Tourniquet | 6.Hand drawn wire saw |
| 7.Scissors | 7.Multi-function survival card |
| 8.Plastic tweezer | 8.Multi-function survival bracelet |
| 9.Pin | 9.Lifeline rope |
| 10.Cotton swab | 10.Multi-function folding knife |



HF CZ01 Compressed Gauze

Product Feature	1. An ultra-compact roll of sterile, high quality 100% cotton gauze, 2. That utilizes a unique crinkle weave that has excellent fluid absorption, 3. Indicated for the treatment of bleeding and the dressing and bandaging of wounds.
Material	100% Cotton
Packing	Sterile, vacuum
Size	4.5"x4.1yards, 6ply



HF BS01 Bandage Scissors

Product Feature	Designed for use by medical professionals or general public by cutting bandages, clothing and any other materials that need immediate access.
Material	Stainless steel
Color	Black or Metallic color



HF EB-1 Emergency Blanket

Product Feature	It is mainly used in hospital operatingroom/field/road traffic to rescue wounded.
Material	Polyester fiber, aluminum coating
Color	Silver, Gold or Green



HF CS-1 Chest Seal

Product Feature	Block inhaled air, release trapped air via valve to prevent tension pneumothorax.
Material	Hydrogel
Type	Vented or Non-vented



HF CAT-1 First Aid Tourniquet

Product Feature	It is mainly used in hospital operatingroom/field/road traffic to rescue wounded.
Color	Black
Size	38x950mm



HF TOU-1 TPE Tourniquet

Product Feature	For the vascular hemorrhage on wounded limbs, the battle line emergency firstaid, also used in pre-hospital emergency treatment
Material	PE
Color	Green; blue etc
Size	1"/2"/3" ect x 0.8/1/1.8/2m etc



HF TB01 Triangular Bandage

Product Feature	1.Immobilizing broken or injured limbs 2.Stopping a hemorrhage by using as a tourniquet with a bar/stick
Material	100% Cotton,100% Non woven
Color	White, Gray or Natural
Size	90x90x127cm,100x100x142cm

急救系列 (FIRST AID)

军用急救包 FIRST AID KIT



HFFAK-E01

- | | |
|-------------------------|-------------------------|
| 1. Compressd Gauze *1 | 7. Chest Seal *1 |
| 2. Emergency Blanket *1 | 8. Emergency Bandage *1 |
| 3. CAT Tourniquet *1 | 9. Bandage Scissors *1 |
| 4. Marker *1 | |
| 5. Alcohol Pad *2 | |
| 6. Elastic Bandage *1 | |



HFFAK-E02

- | | |
|--------------------------|--------------------------|
| 1.Roll Splint *1 | 7.Burn Cream *1 |
| 2.Emergency Blanket *1 | 8.Medical Tape *1 |
| 3.CAT Tourniquet *1 | 9.Elastic Bandage *1 |
| 4.Nasopharyngeal Tube *1 | 10.Chest Seal *1 |
| 5.Marker *1 | 11.Chest Seal (vented)*1 |
| 6.Alcohol Pad *2 | 12.Emergency Bandage *1 |



HFFAK-E03

- | | |
|-----------------------------|-------------------------|
| 1.First Aid Instructions *1 | 7.Medical Tape *1 |
| 2.Compressed Gauze *1 | 8.Elastic Bandage *2 |
| 3.First Aid Blanket *1 | 9.Chest Seal *1 |
| 4.CAT Tourniquet *1 | 10.Emergency Bandage *1 |
| 5.Marker *1 | 11.Bandage Scissors *1 |
| 6.Nasopharyngeal Tube *1 | |



HFFAK-E04

- | | |
|--------------------------|-------------------------|
| 1.Compressed Gauze *1 | 7.Triangular Bandage *1 |
| 2.CAT Tourniquet *1 | 8.Bandage Scissors *1 |
| 3.Marker *1 | |
| 4.Chest Seal *1 | |
| 5.Chest Seal (vented) *1 | |
| 6.Elastic Bandage *1 | |

家用急救包 FIRST HOUSEHOLD KIT



HFFAK-H01

- | | |
|-------------------------|----------------------------|
| 1.PBT Bandage 7.5cm *1 | 6. Non-Adhesive Bandage *2 |
| 2. Medical Tape *2 | 7. Scissors *1 |
| 3. Alcohol Pad *6 | 8. Plastic Tweezer *1 |
| 4. Sterile Gauze Pad *2 | 9. Cotton Swab *15 |
| 5. Adhesive Bandage *3 | |



HFFAK-H02

- | | | |
|-----------------------------|------------------------|------------------------|
| 1. PBT Bandage 5cm *2 | 6. Adhesive Bandage *2 | 11. PVC Gloves *2 |
| 2. PBT Bandage 7.5cm *2 | 7. Non-Adherent Pad *2 | 12. Pin *5 |
| 3. Elastic Bandage 7.5cm *2 | 8. Medical Tape *1 | 13. Scissors *1 |
| 4. Triangular Bandage *1 | 9. Sterile Abd Pad *1 | 14. Plastic Tweezer *1 |
| 5. Alcohol Pad *3 | 10. Cold Pack *2 | |



HFFAK-H03

- | | |
|-------------------------|------------------------|
| 1. PBT Bandage 5cm *1 | 6. Non-Adherent Pad *2 |
| 2. PBT Bandage 7.5cm *1 | 7. Medical Tape *1 |
| 3. Antiseptic Swab *4 | 8. Cold Pack *1 |
| 4. Alcohol Pad *2 | 9. Cotton Swab *10 |
| 5. Adhesive Bandage *3 | |



HFFAK-H04

- | | |
|--------------------------|--------------------|
| 1. PBT Bandage 5cm *1 | 6. Medical Tape *1 |
| 2. Gauze Roll *1 | 7. Scissors *1 |
| 3. Triangular Bandage *1 | 8. Cotton Swab *10 |
| 4. Antiseptic Swab *2 | |
| 5. Adhesive Bandage *5 | |