

PRODUCT DISPLAY



We have a long time experience of foreign trade research and development, we can give the most professional service about the system of Supply From **Demand - Product Design - Mass Production - Transportation Guarantee.**

RUIAN SUNNYWORLD INTERNATIONAL TRADE CO.,LTD. specialize in the production of Tourniquet, Diagnostic Hammer, Flashlight, Stethoscope, Blood Pressure Monitor and other diagnostic products and medical gifts. Our products got CE, ISO, FDA certified and registered with DUNS number and UDI codes. Give customers the most authoritative assurance. In addition, we have product designers, CAD drawings, 3D printed samples. We have three 250T injection molding machines, four 200T injection molding machines, six 150T injection molding machines, and one vertical injection molding machine, which can basically meet the needs of various products. We also have a mold room that can process any mold. Any product can be customized for customer.

At present, we have customized scissors with flashlight, special shape stethoscopes, tourniquets, flashlights, name tags, diagnostic hammers, stethoscope diaphragms, syringe and suture trainers, and other products for customers.

PRODUCT

C/A/T/A/L/O/G

- 01 STETHOSCOPE01
- 12 ANEROID SPHYGMOMANOMETER
- 19 DIGITAL BLOOD PRESSURE MONITOR
- 23 DIGITAL THERMOMETER
- 27 ENT
- 48 TOURNIQUET
- 49 PENLIGHT
- 51 NEUROLOGICAL HAMMER
- 56 MIX SET
- 58 MONOFILAMENT
- 59 OXIMETER
- 63 MEDICAL GIFTS

STETHOSCOPE

SW-ST01A



Single Head Stethoscope

Sliver head(Oxidized&Anodized head for choice)
Standard Y-tubing

For Adult: Diameter of chest piece:46mm

For Child: Diameter of chest piece:37mm

1pc/colour box,100pcs/ctn,52x41x28cm,14/13kgs

SW-ST01B



Colored single head stethoscope

Attractive colored head with matching
color tubing

For Adult: Diameter of chest piece:46mm

For Child: Diameter of chest piece:37mm

1pc/colour box,100pcs/ctn,52x41x28cm,14/13kgs

SW-ST01C



Single Head Stethoscope with non-chill ring

Aluminum chest piece with a wide range of tubing
colors that match non-chill rings

Diameter of chest piece:46mm

For Adult: Diameter of chest piece:46mm

For Child: Diameter of chest piece:37mm

1pc/colour box,100pcs/ctn,52x41x28cm,14/13kgs

SW-ST01F



Plastic Stethoscope

Full set are all produced by plastic

1pc/pp bag,200pcs/ctn,52x41x28cm,13/12kgs

SW-ST02A



Dual head stethoscope

Sliver head(Oxidized&Anodized head for choice)
Standard Y-tubing

For Adult: Diameter of chest piece:46mm

For Child: Diameter of chest piece:37mm

1pc/colour box,100pcs/ctn,52x41x28cm,14/13kgs

SW-ST02B



Colored dual head stethoscope

Attractive colored head with matching color tubing

For Adult: Diameter of chest piece:46mm

For Child: Diameter of chest piece:37mm

1pc/colour box,100pcs/ctn,52x41x28cm,14/13kgs

SW-ST02C



Dual Head Stethoscope with non-chill ring

Aluminum chest piece with a wide range of tubing
colors that match non-chill rings

For Adult: Diameter of chest piece:46mm

For Child: Diameter of chest piece:37mm

1pc/colour box,100pcs/ctn,52x41x28cm,14/13kgs

SW-ST05A



Stethoscope for baby and neonatal use

Chest piece with a wide range of tubing, colors
that match non-chill rings

SW-ST05A1: Aluminum chest piece

SW-ST05A2: Zinc chest piece

SW-ST05A3: Stainless steel chest piece

1pc/colour box,100pcs/ctn,53x40.5x31cm,12/11kgs

STETHOSCOPE

STETHOSCOPE

SW-ST03A



Sprague rappaport stethoscope

Silver color head, standard tube, comes complete with a clear multi-pocket bag for relevant. Accessories and replacement parts. Various color available.

1pc/colour box,50pcs/ctn,66x31x43cm,17/16kgs

SW-ST03B



Sprague rappaport stethoscope

Silver color head, transparent tube, comes complete with a clear multi-pocket bag for relevant

Accessories and replacement parts. Various color available.

1pc/colour box,50pcs/ctn,66x31x43cm,17/16kgs

SW-ST03C



Sprague rappaport stethoscope

Attractive gold-plate head and binaural, more deluxe

comes complete with a clear multi-pocket bag for relevant Accessories and replacement parts. Various color available.

1pc/colour box,50pcs/ctn,66x31x43cm,17/16kgs

SW-ST03D



Sprague rappaport stethoscope

Attractive color head with matching color tubing, comes complete with a clear multi-pocket bag for relevant

Accessories and replacement parts. Various color available.

1pc/colour box,50pcs/ctn,66x31x43cm,17/16kgs

STETHOSCOPE

SW-ST04A



Deluxe Professional Single Head Stethoscope

Zinc alloy chest piece with nickel-plated silver finish.

Comes with special spring binaural tubing.

Special design for get better effect through press the back of head.

Special Design of binaural for suiting the ear structure.

1pc/colour box,50pcs/ctn,61.5x26.5x45cm,15/14kgs

SW-ST04B



Deluxe Professional Single Head Stethoscope

Zinc alloy chest piece with chrome-plated silver finish.

Comes with special spring binaural tubing.

Special design for get better effect through press the back of head.

Special design of binaural for suiting the ear structure.

1pc/colour box,50pcs/ctn,61.5x26.5x45cm,15/14kgs

SW-ST04C



Professional Heart Shape Stethoscope

The professional heart shape stethoscope

made in zinc alloy

heart shape head

carving logo on head or printing logo on diaphragm.

1pc/colour box,50pcs/ctn,61.5x26.5x45cm,15/14kgs

SW-ST04D



Professional Doctor's Single Head Stethoscope

1. Silver single head chest piece , Ergonomic design for the back.so as to get comfortable use. can catch both low frequency and high frequency sound by pressing the back in different pressure.

2. considerate anti-chill ring design

3. connecting spring is embedded ,ear tip is specially designed for suiting the ear structure.

4. zinc alloy chest piece with MATTE SATIN finish.excellect sound conduction. little acoustic attenuation.

1pc/colour box,50pcs/ctn,61.5x26.5x45cm,15/14kgs

STETHOSCOPE

SW-ST06



Special Single Head Stethoscope

Special Single Head Stethoscope.
Zinc head, powder-coating finish.
Comes with non-chill ring.
With inner spring binaural tubing assembly,
Uniquely design of binaural for suiting the ear structure.

1pc/colour box,50pcs/ctn,63x27x43cm,15/14kgs

SW-ST09A



Stethoscope with clock

Aluminum chest piece with a wide range of tubing colors that match non-chill rings,Diameter of chest piece is 47mm.
Comes with inner-spring binaural tubing.
Special design of binaural for suiting the ear structure, the head with clock

1pc/colour box,20pcs/ctn,68x35x25cm,6/5kgs

SW-ST09B



Sprague Rapport Stethoscope with Clock. silver color head, standard tube, with matching color tube, comes complete with a clear multi-purpose bag for relevant. Accessories and replacement parts.Various color a available.

1pc/colour box,20pcs/ctn,68x35x25cm,12/11kgs

SW-ST10A



Sprague rappaport stethoscope for teaching use, with 2 binaural

1pc/colour box,50pcs/ctn,66x31x43cm,25/24kgs

STETHOSCOPE

SW-ST10B



Dual head stethoscope for teaching use, with 2 binaural

1pc/PVC bag,50pcs/ctn,66x34x24cm,15/14kgs

SW-ST10C



Single head stethoscope for teaching use, with 2 binaural

1pc/PVC bag,50pcs/ctn,66x34x24cm,14/13kgs

SW-ST12



Fetal stethoscope for pregnancy use and normal adult use

1pc/PVC bag,50pcs/ctn,66x34x24cm,14/13kgs

SW-ST13



Acrylic stethoscope.

Chest piece in acrylic material with polishing finish inner-spring binaural tubing

1pc/colour box,20pcs/ctn,66x34x24cm,12/11kgs

STETHOSCOPE

SW-ST13B



Nurse heart stethoscope

- 1.copper binaural
- 2.with soft ear tip
- 3.with heart shape red head

1pc/colour box,20pcs/ctn,66x34x24cm,10/9kgs

SW-ST14A



Cardiology master stethoscope

Deluxe zinc single head, special design of head for better sound correct.
 Light copper binaural.
 Chest piece comes with non-chill ring.
 Soft ear-tips with inner spring binaural tubing assembly.
 Special design of binaural for suiting the ear structure, with uniquely designed binaral thick tubing assembly.

1pc/colour box,20pcs/ctn,66x34x24cm,12/11kgs

SW-ST14B



Professional Cardiology Master Stethoscope.

Deluxe zinc single head Special design for better sound correct. Light copper binaural.
 Chest piece comes with non-chill ring. Soft ear-tips with inner spring binaural tubing assembly. Special design of binaural for suiting the ear structure, with uniquely designed binaural thick tubing assembly.

1pc/colour box,20pcs/ctn,66x34x24cm,12/11kgs

SW-ST14E



Professional Cardiology Master Stethoscope. Heart Shape

- 1.Heart shape chest piece , Ergonomic design for the back.so as to get comfortable use.
can catch both low frequency and high frequency sound by pressing the back in different pressure.
- 2.with integrated diaphragm (ruber ring and diaphragm is one set)
- 3.connecting spring is embedded ,ear tip is specially designed for suiting the ear structure.
- 4.zinc alloy chest piece with Nickel-plated or Chrome-plated finish.excellect sound conduction. little acoustic attenuation .

1pc/colour box,20pcs/ctn,66x34x24cm,12/11kgs

STETHOSCOPE

SW-ST15A



Deluxe stainless steel head,on side for child, the other side for adult.

Light copper binaural.

Chest piece comes with non-chill ring, soft ear-tips with inner spring binaural tubing assembly.

Special design of binaural for suiting the ear structure,with uniquely designed binaural thick tubing assembly.

1pc/colour box,20pcs/ctn,66x34x24cm,12/11kgs

SW-ST16



Deluxe stainless steel stethoscope

Light copper binaural. Chest piece comes with non-chill ring, soft ear-tips, with inner spring binaural tubing assembly. Special design of binaural for suiting the ear structure.

For Adult: Diameter of chest piece:46mm

For Child: Diameter of chest piece:37mm

1pc/colour box,20pcs/ctn,66x34x24cm,12/11kgs

SW-ST18



Cardiology stainless steel stethoscope,for adult,cardiology type,feature a deep cone-shaped bell.
 Chest piece comes with non-chill ring.

With inner spring binaural tubing assembly.

Special design of binaural for suiting the ear structure.

With uniquely designed binaural thick tubing assembly.

1pc/colour box,20pcs/ctn,66x34x24cm,12/11kgs

SW-ST19



Cardiology stainless steel stethoscope, Stainless steel stethoscope.

Three Head (For adult,For child,For baby),

Triangle-shaped tube with Stainless Steel Binaural.

1pc/colour box,20pcs/ctn,71x39x25cm,14.5/13.5kgs

STETHOSCOPE

SW-ST21



Stainless steel single head stethoscope, for adult.
Deluxe single head style.
With inner spring binaural tubing assembly.
Special design of binaural for suiting the ear structure.

1pc/colour box, 20pcs/ctn, 66x34x24cm, 12/11kgs

SW-ST24



Deluxe dual head stethoscope.
Features a developed new design chrome-plated zinc alloy chest piece screw-on retaining ring, and internal spring binaural, matched tubing, and non-chill ring, for professional purpose.

1pc/colour box, 50pcs/ctn, 70x43x28cm, 17/16kgs

SW-ST27A



Cardiology stethoscope

Zinc alloy chest piece, copper, PVC, coated by Chrome

1pc/colour box, 20pcs/ctn, 66x34x24cm, 12/11kgs

SW-ST27B



Cardiology stethoscope

Zinc alloy chest piece

1pc/colour box, 20pcs/ctn, 66x34x24cm, 12/11kgs

STETHOSCOPE

SW-ST27C



Cardiology stethoscope

With black color plating, zinc alloy chest piece

1pc/colour box, 20pcs/ctn, 66x34x24cm, 12/11kgs

SW-ST33



Master cardiology type stethoscope

Zinc alloy chest piece

1pc/colour box, 20pcs/ctn, 66x34x24cm, 12/11kgs

SW-ST47




Cardiology stainless steel stethoscope

With extra different shape chest piece, can be converted
Gift foam box packing or carry case packing

1pc/colour box, 25pcs/ctn, 84x33x24cm, 16/15kgs

SW-ST48



Cardiology stainless steel stethoscope

The height: 1.9cm
Material: Stainless steel
Color: Dazzling colors color, golden color, Silver color, Rose gold color
Diaphragm and rubber ring as a whole, frequency conversion
Ear tips: silicon
Tube: can offer any color and transparent
With extra accessory bag: 2pcs ear tips, 1pc diaphragm
Packing: gift foam box

1pc/colour box, 25pcs/ctn, 84x33x24cm, 16/15kgs

STETHOSCOPE

SW-AP22



Menbrane and Ruber as a whole ,Any color viable

SW-AP23



Carry case for stethoscope and other items, any color viable

SW-G03A



Wall type stethoscope holder. Paste the holder on the wall

SW-G03B



Hip-clip stethoscope holder. Hang the holder on the waist

SW-G02Y




Cover shape stethoscope ID tag

SW-G02A



Stethoscope ID tag

SW-G02N



Stethoscope ID tag

SW-G02E



Stethoscope ID tag

SW-G02R



Stethoscope ID tag

SW-02Q



Stethoscope ID tag

SW-AP02



Binaural

SW-G02R



Stethoscope tubing

ANEROID SPHYGMOMANOMETER

SW-AS01



Cuff: calibrated nylon or cotton cuff without metal ring
Gauge: stop pin or non-stop pin, 120g or 145g
Material of gauge: aluminium or zinc alloy
standard valve and bulb
vinly zipper case or nylon bag packing
SW-AS01A: with latex bulb and latex bladder
SW-AS01B: with PVC bulb and PVC bladder

Packing: 1pc/ zipper bag/ box,50pcs/ctn,
56x39x40cm,20/19kgs

SW-AS02



Cuff: standard nylon or cotton with metal-ring cuff
Gauge: stop pin or non-stop pin, 120g or 145g
Material of gauge: aluminium or zinc
standard valve and bulb
vinly zipper case or nylon bag packing
SW-AS02A: with latex bulb and latex bladder
SW-AS02B: with PVC bulb and PVC bladder

Packing: 1pc/ zipper bag/ box,50pcs/ctn,
56x39x40cm,20/19kgs

SW-AS03



Cuff: nylon or cotton cuff with metal hook
Gauge: stop pin or non-stop pin, 120g or 145g
Material of gauge: aluminium or zinc
standard valve and bulb
vinly zipper case or nylon bag packing
SW-AS03A: with latex bulb and latex bladder
SW-AS03B: with PVC bulb and PVC bladder

Packing: 1pc/ zipper bag/ box,50pcs/ctn,
56x39x40cm,20/19kgs

SW-AS04



Cuff: nylon or cotton cuff with metal chip vinly
zipper case or nylon bag packing
Gauge: stop pin or non-stop pin, 120g or 145g
Material of gauge: aluminium or zinc
SW-AS04A: with latex bulb and latex bladder
SW-AS04B: with PVC bulb and PVC bladder

Packing: 1pc/ zipper bag/ box,50pcs/ctn,
56x39x40cm,20/19kgs

ANEROID SPHYGMOMANOMETER

SW-AS05



Cuff: D-ring cuff or calibrated cuff
 Gauge: stop pin or non-stop pin,
 Material of gauge: zinc ally + plastic case
 standard valve and bulb
 vinly zipper case or nylon bag packing
 SW-AS05A: with latex bulb and latex bladder
 SW-AS05B: with PVC bulb and PVC bladder

Packing: 1pc/ zipper bag/ box,50pcs/ctn,
 56x39x40cm,20/19kgs

SW-AS06



nylon or cotton cuff with leather attachment
 Gauge: Fasten to the cuff, stop pin or non-stop pin,
 chrome plated manometer with rubber
 protective ring.
 SW-AS06A: with latex bulb and latex bladder
 SW-AS06B: with PVC bulb and PVC bladder

Packing: 1pc/ zipper bag/ box,50pcs/ctn,
 56x39x40cm,20/19kgs

SW-AS07



Palm type sphygmomanometer
 Gauge: chrome plated plastic case , with single
 tube or with double tube
 SW-AS07A: with latex bulb and latex bladder
 SW-AS07B: with PVC bulb and PVC bladder

Packing: 1pc/zipper bag/colour box,50pcs/ctn,
 56x39x40cm,25/24kgs

SW-AS08



Palm type sphygmomanometer
 Gauge: with single tube or with double tube, ABS
 plastic case, stop pin or non-stop pin,
 black powder-coating large handle
 SW-AS08A: with latex bulb and latex bladder
 SW-AS08B: with PVC bulb and PVC bladder

Packing: 1pc/zipper bag/colour box,50pcs/ctn,
 56x39x40cm,25/24kgs

ANEROID SPHYGMOMANOMETER

SW-AS09



Palm type sphygmomanometer
 Gauge: with single tube or with double tube, ABS
 plastic case, stop pin or non-stop pin,black
 powder-coating large handle
 SW-AS09A: with latex bulb and latex bladder
 SW-AS09B: with PVC bulb and PVC bladder

Packing: 1pc/zipper bag/colour box,50pcs/ctn,
 56x39x40cm,25/24kgs

SW-AS10



Palm type sphygmomanometer
 Gauge: with single tube or with double tube, chrome
 plated plastic case, 50mm diameter dial.
 Stainless steel handle
 SW-AS10A: with latex bulb and latex bladder
 SW-AS10B: with PVC bulb and PVC bladder

Packing: 1pc/zipper bag/colour box,50pcs/ctn,
 56x39x40cm,25/24kgs

SW-AS11



Palm type sphygmomanometer
 Gauge: Zinc alloy with black coating, 65mm
 diameter dial. With single tube or with
 double tube
 SW-AS11A: with latex bulb and latex bladder
 SW-AS11B: with PVC bulb and PVC bladder

Packing: 1pc/zipper bag/colour box,50pcs/ctn,
 56x39x40cm,25/24kgs

SW-AS12



Palm type sphygmomanometer
 Gauge: zinc alloy with black coating, 50mm
 diameter dial, Black stainless steel handle.
 With single tube or with double tube

Packing: 1pc/zipper bag/colour box,50pcs/ctn,
 56x39x40cm,25/24kgs

ANEROID SPHYGMOMANOMETER

SW-AS13



Aneroid sphygmomanometer with Single head stethoscope

Packing: 1set/zipper bag/colour box,50sets/ctn, 56x39x40cm,25/24kg

SW-AS14



Aneroid sphygmomanometer with Dual head stethoscope

Packing: 1set/zipper bag/colour box,50sets/ctn, 56x39x40cm,25/24kg

SW-AS15



Aneroid sphygmomanometer with rappaport stethoscope

1set/carrying bag/colour box,20sets/ctn, 47.5x29x31.5cm,14/13kgs

SW-AS16



Desk/Wall type sphygmomanometer
Material: ABS plastic
Square shaped 14cm diameter dial front and rear adjustable for easy reading.
with or without PVC coil tubing.

packing: 1pc/colour box,20pcs/ctn,85x39x38cm, 22/21kgs

ANEROID SPHYGMOMANOMETER

SW-AS17



Desk type sphygmomanometer
ABS plastic round dia with metal basket.
With or without PVC coil tubing.

Packing: 1pc/box,6pcs/ctn,68x37x37cm, 8.5/7.5kgs

SW-AS18



DEXLUXE Palm type sphygmomanometer
Gauge: Zinc case with chrome-plated , 65mm diameter dial.
PVC ribbed bulb.

Packing: : 1pc/zipper bag/colour box,50pcs/ctn, 56x39x40cm,25/24kgs

SW-AS19



Standing type sphygmomanometer
Gauge: Square shaped 14cm diameter dial,360 manometer swivel for easy-reading
Front and rear adjustable
5-wheel with brake with aluminium rod, adjustable high(with range of 85-125cm)
with or without PVC coil tubing

Packing: 1pc/ctn,80x32x17cm,4/3kgs

SW-AS20



Cuff: nylon or cotton cuff with metal chip
vinyl zipper case or nylon bag packing
Gauge: stop pin or non-stop pin, 120g or 145g
Material of gauge: aluminium or zinc
SW-AS04A: with latex bulb and latex bladder
SW-AS04B: with PVC bulb and PVC bladder

Packing: 1pc/ zipper bag/ box,50pcs/ctn, 56x39x40cm,20/19kgs

ANEROID SPHYGMOMANOMETER

SW-AS21



Standing type blood pressure monitor
360 sphygmomanometer switch
5-wheel with brake, Aluminum tube, adjustable high stand
(With range of 95-135cm)
Nylon cuff with D-ring
2-tube Adult Latex Bulb
Standard Latex Bulb
Air-release Valve with Spring
Standard End Valve
PVC coil tubing

Packing: 1pc/CTN, 80x32x17cm, 5/3kgs

SW-AS24



Sphygmomanometer with stethoscope
The stethoscope chest piece is screwed to cuff

Packing: 1set/zipper bag/colour box, 50sets/ctn,
56x39x40cm, 25/24kgs

SW-AS25



Palm type sphygmomanometer
Plastic ABS case
64.5mm-diameter dial
Three mode of handle for choice One-whole-set
deflation valve
Nylon cuff with D-ring
1-tube PVC bladder
Flat bottom PVC bulb
Taper metal net end valve

Packing: 1pc/carrying bag/colour box, 30pcs/ctn,
48x44x32.5cm, 17/16kgs

SW-AS27



Palm type sphygmomanometer
Classic shining nylon cuff without D-ring for child/
adult/obese adult
1-tube bladder for child/ adult/obese adult
Standard latex bulb
Standard end valve
Shining nylon carrying bag

Packing: : 1set/ box, 20pcs/ctn, 56x50x30.5cm,
17.5/16kgs

ANEROID SPHYGMOMANOMETER

SW-AB01



SW-AB02A



SW-AB02B



SW-AB02C



SW-AB03A



SW-AB03B



SW-AB04A



SW-AB04B



SW-AB04C



SW-AB05A



SW-AB05B



SW-AB06



SW-AB08



DIGITAL BLOOD PRESSURE MONITOR

SW-DBP80E



Arm Blood Pressure Monitor
 LCD size: 55×44mm with digit height 18mm
 Product size: 132×100×45mm
 • WHO indicator • 2×90 memory
 • 3 times average • IHB indicator
 • Year/month/date/time function
 • Beep alarm • Automatic power -off
 • Power by 4 AA size alkaline batteries or AC/DC power supply

SW-DBP60BH



Wrist Blood Pressure Monitor
 LCD size: 47×39mm with digit height 14mm
 Product size: 79×66×64.5mm
 • IHB indicator
 • WHO indicator
 • 2×90 memory
 • 3 times average
 • Year/month/date/time function
 • Automatic power off
 • Power by 2 AAA size alkaline batteries

SW-DBP60CH



Wrist Blood Pressure Monitor
 LCD size: 47×39mm with digit height 14mm
 Product size: 79×66×64.5mm
 • IHB indicator
 • WHO indicator
 • 2×90 memory
 • 3 times average
 • Year/month/date/time function
 • Automatic power off
 • power by 2 AAA size alkaline batteries

SW-DBP80EH



Arm Blood Pressure Monitor
 LCD size: 76.8×60mm with digit height 22mm
 Product size: 132×100×45mm
 • WHO indicator • 2×90 memory
 • 3 times average • IHB indicator
 • Year/month/date/time function
 • Beep alarm • Automatic power off
 • Power by 4 AA size alkaline batteries or AC/DC power supply

DIGITAL BLOOD PRESSURE MONITOR

SW-DBP80K



Arm Blood Pressure Monitor
 LCD size: 80×72mm
 Product size: 139×100×46mm
 • Blue backlight for LCD & 5START/STop button
 • WHO indicator • 2×90 memory
 • 3 times average • Year/month/date/time function
 • IHB indicator. Beep alarm. Automatic power-off
 • power by 4 AA size alkaline batteries or AC/DC power supply

SW-DBP802



Arm Blood Pressure Monitor
 • Smart 5sense
 • WHO Indication
 • Power Port(optional)
 • 2x120 MemfIPs(2 Users)
 • Heart Pulse Arrhythmia Detection
 • Average For Last 3 Measurements

SW-DBP807-1



Arm Blood Pressure Monitor
 • 5matt Sense
 • WHO Indication
 • Micro USB (optional)
 • 2x120 Memories(2 Users)
 • Heart Pulse Arrhythmia Detection
 • Average For Last 3 Measurements

SW-DBP811



Arm Blood Pressure Monitor
 • Smart Sense
 • WHO Indication
 • Power Port (optional)
 • 2x120 Memories(2 Users)
 • Heart Pulse Arrhythmia Detection
 • Average For last 3 Measurements

DIGITAL BLOOD PRESSURE MONITOR

SW-DBP806-1



Arm Blood Pressure Monitor

- Bluetooth
- Smart Sense
- WHO Indication
- Micro USB (optional)
- Body motion detection
- 2x120 Memories(2 Users)
- Heart Pulse Arrhythmia Detection
- Average For Last 3 Measurements
- Automatic detection of cuff tightness

SW-DBP801



Arm Blood Pressure Monitor

- Smart Sense
- WHO Indication
- Power Port (optional)
- 2x120 Memories(2 Users)
- Heart Pulse Arrhythmia Detection
- Average For last 3 Measurements

SW-DBP-805



Arm Blood Pressure Monitor

- Smart Sense
- WHO Indication
- Power Port (optional)
- 2x120 Memories(2 Users)
- Heart Pulse Arrhythmia Detection
- Average For last 3 Measurements

SW-DBP3

Arm blood pressure monitor



Arm Blood Pressure Monitor

- Unit conversion: mmHg, kPa/mmHg
- 2 users, 99 memories per user
- IHB, Irregular heart beating sign
- WHO, indicating blood pressure level in 4 colors
- Average of latest 3 records
- Low power detection
- Over pressure protection, 295mmHg (20ms)
- Automatic switch off
- Error message indication, date & time indication
- Voice (option), AC adaptor (option)

DIGITAL BLOOD PRESSURE MONITOR

SW-DBP4

Arm blood pressure monitor



Arm Blood Pressure Monitor

- Large LCD display: 72mm x 43.5mm
- Unit conversion: mmHg, kPa/mmHg
- 2 users, 99 memories per user
- IHB, Irregular heart beating sign
- WHO, indicating blood pressure level in 4 colors
- Average of latest 3 records
- Low power detection
- Over pressure protection, 295mmHg (20ms)
- Automatic switch off
- Error message indication, date & time indication
- Voice (option), AC adaptor (option)

SW-DBPW2

Wrist blood pressure monitor



Wrist Blood Pressure Monitor

- Large LCD display: 33.8mm x 29mm
- Unit conversion: mmHg, kPa/mmHg
- 2 users, 99 memories per user
- IHB, Irregular heart beating sign
- WHO, indicating blood pressure level in 4 colors
- Average of latest 3 records
- Low power detection
- Over pressure protection, 295mmHg (20ms)
- Automatic switch off
- Error message indication, date & time indication
- Voice (option), AC adaptor (option)

DIGITAL THERMOMETER

SW-DT01



Measured temperature range: 32 centigrade ---42 centigrade
 Centigrade will display when the temperature is below 32 centigrade
 Centigrade will display when the temperature is above centigrade
 Measurement precision: ± 0.1 centigrade
 (35 centigrade - 39 centigrade) ± 2 centigrade
 (below 35 centigrade and above 39 centigrade)
 Display mode: LCD (three and a half bits)
 Power supply mode: AG3 type 1.5V normal button battery
 User hint: After completion of the measurement,
 buzz sound will be given out.

SW-DT01 B



Measured temperature range: 32 centigrade ---42 centigrade
 Centigrade will display when the temperature is below 32 centigrade
 Centigrade will display when the temperature is above centigrade
 Measurement precision: ± 0.1 centigrade
 (35 centigrade - 39 centigrade) ± 2 centigrade
 (below 35 centigrade and above 39 centigrade)
 Display mode: LCD (three and a half bits)
 Power supply mode: AG3 type 1.5V normal button battery
 User hint: After completion of the measurement,
 buzz sound will be given out.

SW-DT02



Measured temperature range: 32 centigrade ---42 centigrade
 Centigrade will display when the temperature is below 32 centigrade
 Centigrade will display when the temperature is above centigrade
 Measurement precision: ± 0.1 centigrade
 (35 centigrade - 39 centigrade) ± 2 centigrade
 (below 35 centigrade and above 39 centigrade)
 Display mode: LCD (three and a half bits)
 Power supply mode: AG3 type 1.5V normal button battery
 User hint: After completion of the measurement,
 buzz sound will be given out.

SW-DT03



Range: Centigrade 32.0-42.0°C
 Fahrenheit 90.0-108.0°F
 Accuracy : Centigrade $\pm 0.2^\circ\text{C}$ less than 35.0°C
 $\pm 0.1^\circ\text{C}$ 35.0°C to 39.0°C
 $\pm 0.2^\circ\text{C}$ Greater than 39.0°C
 Fahrenheit $\pm 0.3^\circ\text{F}$ less than 95.0°F
 $\pm 0.2^\circ\text{F}$ 95.0°F to 102.0°F
 $\pm 0.3^\circ\text{F}$ Greater than 102.0°F
 At standard room temperature 25.0°C (77.0°F)
 Display: Liquid crystal display 3 1/2 digits
 Battery: One 1.55V button size battery (UCC 392 or LR41)
 Power consumption: 0.15mW in measurement mode
 Battery life: More than 200 hours of continuous operation
 Weight: Approx 12 grams including battery
 Dimensions: 12mm by 22mm by 145mm
 Alarm (beeper mode only) Approx 10 second signal when peak temperature reached

DIGITAL THERMOMETER

SW-OT04



- 1) Display range: 32.0-42.0°C (90-107.6°F)
- 2) Accuracy: $\pm 0.1^\circ\text{C}$ ($\pm 0.2^\circ\text{F}$)
- 3) Min. scale: 0.1
- 4) Measurement time (reference only): a) 60 $\pm$ 10 seconds (oral) b) 100 $\pm$ 20 seconds (underarm)
- 5) Beeper function
- 6) Auto shut-off
- 7) Battery: 1.5V button battery (LR/SR-41)
- 8) Size: 126 x 18 x 10mm
- 9) LCD: 20 x 7mm
- 10) N.W.: 11.5g
- 11) waterproof (or not)

SW-DT04B



- 1) Display range: 32.0-42.0°C (90-107.6°F)
- 2) Accuracy: $\pm 0.1^\circ\text{C}$ ($\pm 0.2^\circ\text{F}$)
- 3) Min. scale: 0.1
- 4) Measurement time (reference only):
a) 60 $\pm$ 10 seconds (oral) b) 100 $\pm$ 20 seconds (underarm)
- 5) Beeper function
- 6) Auto shut-off
- 7) Battery: 1.5V button battery (LR/SR-41)
- 8) Size: 126 x 18 x 10mm
- 9) LCD: 20 x 7mm
- 10) N.W.: 11.5g
- 11) waterproof (or not)

SW-DT04C



- 1) Display range: 32.0-42.0°C (90-107.6°F)
- 2) Accuracy: $\pm 0.1^\circ\text{C}$ ($\pm 0.2^\circ\text{F}$)
- 3) Min. scale: 0.1
- 4) Measurement time (reference only):
a) 60 $\pm$ 10 seconds (oral) b) 100 $\pm$ 20 seconds (underarm)
- 5) Beeper function
- 6) Auto shut-off
- 7) Battery: 1.5V button battery (LR/SR-41)
- 8) Size: 126 x 18 x 10mm
- 9) LCD: 20 x 7mm
- 10) N.W.: 11.5g
- 11) waterproof (or not)

SW-DT10A,B,C,D,E,F,G



Flexible Tip cartoon Type digital thermometer
 Measuring range : 32.0°C - 43.0°C Fahrenheit and
 Centigrade switchable Accuracy: $\pm 0.1^\circ\text{C}$
 Unit size: L 125MM X W 19MM X H 9MM Weight: 9.5g
 Safe, fast and professional accuracy
 Delivers reading in 60 seconds
 3-in-1 Functions: Takes oral, rectal,
 and underarm temperatures
 LCD display; easy reading
 Memory recall for last reading
 Automatic shut-off to help
 extend battery life
 Include storage case
 Certificate: ISO 13485, RoHS, CE 1023

DIGITAL THERMOMETER

SW-DT06A



Nipple Thermometer for baby use
packaging: 1pc/box

SW-DT16



Excellent adaption to ambient temperature.
Accurat reliable even udner complicated surrounding.
Body mode and object mode available.
Warning for fever or high temperature and 3 color back light(Red, Yellow, Green, alarm value can defined by user)
CjF reading available
Auto shut-off: 10 seconds without operation

SW-DT121



Features
Excellent adaption to ambient temperature.
Accurat reliable even udner complicated surrounding.
Body mode and object mode available.
Warning for fever or high temperature and 3 color back light(Red, Green, alarm value can defined by user)
CjF reading available
Auto shut-off: 10 seconds without operation

SW-DT15A2



Features
Excellent adaption to ambient temperature
Accurat reliable even udner complicated surrounding.
Body mode and object mode available. One switch change mode.
Beeper for fever or high temperature (alarm value can defined by user)
C/F reading available
Auto shut-off: 10 seconds without operation

DIGITAL THERMOMETER

SW-DT123



Features
Excellent adaption to ambient temperature.
Accurat reliable even udner complicated surrounding
Body mode and object mode available. One switch change mode.
Beeper for fever or high temperature (alarm value can defined by user)
C/F reading available
Auto shut-off: 10 seconds without operation
No laser point

SW-DT8806c



Features
Excellent adaption to ambient temperature.
Accurat reliable even udner complicated surrounding.
Body mode and 0 均ect mode available. One switch change mode.
Beeper for fever or high temperature (alarm value can defined by user)
C/F reading available
Auto shut-off: 10 seconds without operation

SW-DTH10





Infrared thermometer
Temperature of body body or surface of objecte
Detailed 34 memories with time/date
Response time: 1 second
Effective distance: 1-15cm
User guidance by infrared beam
Fever alarm sounds
Selectable °C/F mode
Error massage indication
Fever teller
150mm x 75mm x 50mm , 120g
2 x AAA batteries

SW-DTH9



The infrared thermometer we produced is specially for measuring body temperature, it can measure human ear and forehead temperature, it suitable for medical unit and home use.
Otitis externa, otitis media
1) This machine has two functions, it can measure both ear and forehead temperature.
2) 1 second measure the temperature, easily and fast.
3) Sensor measurement technology, high precision.
4) Automatically power-off, if left idle for 60 seconds.
5) One-key measurement, easy to use.
6) Alarm for fever, better to know your body situation.
7) Stores 12 sets recent measurement data, easy for your data contrast.
8) Safety by infrared measuring, avoid the damage of the measuring by traditional mercury thermometer.



Specification

- Direct illumination
- 3X magnification for the viewing window
- 2.4mm/3mm/4mm/5mm ear tips of two groups
- The clip with integrated on/off switch
- 3200k and 6000k color temperature LED bulb and XHL Xenon Halogen bulb to optional
- Colors to choose



- 直接照明
- 3倍放大观察窗
- 2.4mm/3mm/4mm/5mm耳套2组
- 固定夹扣式on/off 开关
- 3200k和6000k色温LED灯和XHL卤素灯可选
- 五种颜色可选择



Direct illumination

SW-OT16J

Package content

Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
193x102x32	400x215x195	20	4.9/4.5

Gross and net weight(Kgs/Ctn)Battery containing

Specification

- Direct illumination
- 3X magnification for the viewing window
- 2.4mm/3mm/4mm/5mm ear tips of two groups
- 3200k and 6000k color temperature LED bulb and XHL Xenon Halogen bulb to optional
- Standard handle, handy and convenient
- The clip with integrated on/off switch

- 直接照明
- 3倍放大观察窗
- 2.4mm/3mm/4mm/5mm耳套2组
- 3200k和6000k色温LED灯和XHL卤素灯可选
- 标准型手柄，轻巧方便
- 固定夹扣式on/off 开关



Direct illumination

Package content

Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
190x110x40	400x240x230	20	6.8/6.3

Gross and net weight(Kgs/Ctn)Battery not included



SW-OT16G



SW-OT16N

Specification

- Fiber optic illumination
- 3X magnification for the viewing window
- 2.4mm/3mm/4mm/5mm ear tips of two groups
- 3200k and 6000k color temperature LED bulb and XHL Xenon Halogen bulb to optional
- Supplied by 2xAA alkaline batteries
- 光纤照明
- 3倍放大观察窗
- 2.4mm/3mm/4mm/5mm耳套2组
- 3200k和6000k色温LED灯和XHL卤素灯可选
- 2×AA碱性电池



Fiber optic illumination


Package content

Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
193x102x32	400x215x195	20	4.9/4.5

Gross and net weight(Kgs/Ctn)Battery containing

Specification

- Fiber optic illumination
- 3X magnification for the large size viewing window
- 2.4mm/3mm/4mm/5mm ear tips of two groups
- 3200k and 6000k color temperature LED bulb and XHL Xenon Halogen bulb to optional
- The socket of ear insufflation bulb
- 光纤照明
- 3倍放大的大观察窗
- 2.4mm/3mm/4mm/5mm耳套2组
- 3200k与6000k色温LED灯和XHL卤素灯可选
- 气囊球插口



Fiber optic illumination

Package content

Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
193x102x32	400x215x195	20	4.5/4.1

Gross and net weight(Kgs/Ctn)Battery not included





SW-OT16S



SW-OT16

Specification

- Fiber optic illumination
- 3X magnification for the viewing window
- 2.4mm/3mm/4mm/5mm ear tips of two groups
- 3200k and 6000k color temperature LED bulb and XHL Xenon Halogen bulb to optional
- 光纤照明
- 3倍放大观察视窗
- 2.4mm/3mm/4mm/5mm耳套2组
- 3200k和6000k色温的LED灯和XHL卤素灯可选



Fiber optic illumination

Package content


Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
190x110x40	400x240x230	20	6.8/6.3



Color optional

Specification

- Fiber optic illumination
- 3X magnification for the large size viewing window
- 2.4mm/3mm/4mm/5mm ear tips of two groups
- The socket of ear insufflation bulb
- 3200k and 6000k color temperature of LED bulb and XHL xenon halogen bulb to optional
- 光纤照明
- 3X放大的大观察视窗
- 2.4mm/3mm/4mm/5mm耳套2组
- 气囊球插口
- 3200k和6000k色温的LED灯和XHL卤素灯可选



Fiber optic illumination

Package content

BOX SIZE(MM)	CARTON SIZE(MM)	PCS/CTN	KGS/CTN
190x110x40	400x240x230	20	7.0/6.5

Gross and net weight(Kgs/Ctn)Battery not included




SW-OT16D



SW-OT16L

Specification

- 3X magnification for the large size viewing window
- Fiber optic illumination
- The socket of ear insufflation bulb
- The deluxe LCD screen handle can adjust the brightness of the light source
- Optical SMT LED technology, low power consumption and long life
- CRI≥90%, It's more color reducibility
- Built-instant current IC circuit, The light effect is more stable
- 3X大观察视窗
- 光纤照明系统
- 气囊球插口
- 豪华液晶屏手柄，可以调节光源亮度
- Optical SMT LED技术，低功耗、长寿命
- 显指≥90%色彩还原性强
- 手柄内置恒流稳压芯片，光效稳定



FIBER OPTIC ILLUMINATION

Package content

Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
190x110x40	400x240x230	20	7.0/6.5

Specification

- Anti-scratch glass large viewing window of 3x magnification Fiber optic illumination
- The socket of ear insufflation bulb
- Supplied by 2xAA alkaline batteries (not included)
- The deluxe LCD screen handle can adjust the brightness of the light source
- 防刮玻璃3倍放大的大观察窗
- 光纤照明系统
- 2×AA碱性电池（不包含）
- 豪华液晶屏手柄，可以调节光源亮度



GIFT BOX OF INSUFFLATION BULB



SW-OT16R

Package content

Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
200x140x57	420x300x320	20	7.5/7.0

Gross and net weight(Kgs/Ctn)Battery not included



Package content

Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
200x140x57	420x300x320	20	11.1/10.6

Gross and net weight(Kgs/Ctn)Battery containing



SW-OT16R

SW-OT16F


Specification

- Fiber optic illumination
- Anti-scratch glass large viewing window of 3x magnification
- 3200k or 6000k color temperature LED bulb as light source
- For pneumatic testing
- Li-ion rechargeable battery, long life
- Micro USB charger interface is very convenient for charging

- 光纤照明
- 防刮玻璃3倍放大的大观察窗
- 3200k或6000k色温LED灯光源
- 可鼓膜气动检查
- 锂电池供电，长寿命
- MicroUSB充电接口，充电方便



Two charging methods



Fiber optic illumination



SW-OT15S

Specification


- Simple handle
- 5 different apertures
- 19 kinds of adjustable corrective lenses, Range of lenses : -20d~+20d
- 4000k color temperature LED bulb and XHL Xenon Halogen bulb to optional
- The clip with integrated on/off switch
- 简易型手柄
- 5种光阑
- 19种屈光透镜, 屈光范围: -20D~+20D
- 4000k色温LED灯和XHL卤素灯可选
- 固定夹扣式on/off开关

Package content

Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
193x102x32	400x215x195	20	4.0/3.6

Gross and net weight(Kgs/Ctn)Battery not included

+	1	2	3	4	6	8	10	15	20
-	1	2	3	4	6	8	10	15	20



Specification

- Coaxial parallel light principle
- 5 different apertures
- 19 kinds of adjustable corrective lenses, Range of lenses : -20d~+20d
- 4000k color temperature LED bulb and XHL Xenon Halogen bulb to optional
- The clip with integrated on/off switch
- Standard handle, handy and convenient
- 同轴平行光原理
- 5种光阑
- 19种屈光透镜, 屈光范围: -20D~+20D
- 4000k色温LED灯和XHL卤素灯可选
- 固定夹扣式on/off开关
- 标准型手柄, 轻巧方便



Package content

Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
190x110x40	400x240x230	20	6.5/6.0

Gross and net weight(Kgs/Ctn)Battery not included



SW-OT15




SW-OT15R

Specification

- Optical glass lens and 5 different apertures
- 19 kinds of adjustable corrective lenses, Range of lenses : -20d~+20d
- Supplied by 2xAA alkaline batteries (not included)
- The deluxe LCD screen handle can adjust the brightness of the light source
- Optical SMT LED technology, low power consumption and long life
- Built-instant current IC circuit, The light effect is more stable
- 3200k color temperature of warm white light
- CRI≥90%, It's more color reducibility to ensure that the real image reproduction

- 光学玻璃透镜, 五种不同光阑
- 19种屈光透镜, 屈光范围: -20D~+20D
- 2xAA碱性电池 (不包含)
- 豪华液晶屏手柄, 可以调节光源亮度
- Optical SMT LED技术, 低功耗、长寿命
- 手柄内置恒流稳压芯片, 光效稳定
- 3200K色温的暖白光照明
- 色彩还原性强, 显指≥90%保证图像真实



Package content

Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
190x110x40	400x240x230	20	6.5/6.0

Gross and net weight(Kgs/Ctn)Battery not included

+	1	2	3	4	6	8	10	15	20
-	1	2	3	4	6	8	10	15	20

SW-OT15F

Specification

- Coaxial parallel light principle
- 6 different apertures with red-free filter & blue filter they can be switched
- 29 kinds of adjustable corrective lenses, Range of lenses : -35D~+40D
- Supplied by 2XAA alkaline batteries (not included)
- Deluxe LCD screen handle, can be adjustable the brightness of light source
- CRI≥90%It's more color reducibility
- Constant current IC circuit, maintain the optimal lighting effect
- Low power consumption, long life of 20000 h

- 同轴平行光原理
- 6种光阑和无赤滤光、钴蓝滤光分别与每个光阑切换
- 29种屈光透镜, 屈光范围: -35D~+40D
- 2×A碱性电池 (不包含)
- 豪华液晶屏手柄, 可以调节光源亮度
- 色彩还原性强, 显指≥90%
- 内置恒流稳压芯片, 光效稳定
- 低功耗LED灯, 寿命长达2万小时



Package content

Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
190x100x40	400x240x230	20	6.2/5.7

Gross and net weight(Kgs/Ctn) Battery not included



**Package content**

Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
275x160x56	340x295x310	10	9.8/9.1

Gross and net weight(Kgs/Ctn)Battery containing

- The handle can be selected as a standard type of handle
- 手柄可选配为标准型手柄



**SW-OT17SA**

- Simple handle
- 2.4mm/3mm/4mm/5mm ear tips of one groups
- A head of fiber optic otoscope (OT10S)
- A head of direct ophthalmoscope (OP10S)
- 简易型手柄
- 两组2.2mm/3mm/4mm/5mm 耳套
- 一个光纤检耳镜 (OT10S) 头部
- 一个直接检眼镜 (OP10S) 头部

**SW-OT17**

- A standard handle
- 2.4mm/3mm/4mm/5mm ear tips
- A head of fiber optic otoscope (OT10)
- A head of direct ophthalmoscope (OP10)
- 一个标准型手柄
- 一组2.2mm/3mm/4mm/5mm 耳套
- 一个光纤检耳镜 (OT10) 头部
- 一个直接检眼镜 (OP10) 头部

**SW-OT17C**

- A fiber optic otoscope (OT10)
- A direct ophthalmoscope (OP10)
- 2.4mm/3mm/4mm/5mm ear tips two groups
- 一个光纤检耳镜 (OT10)
- 一个直接检眼镜 (OP10)
- 两组2.2mm/3mm/4mm/5mm 耳套

**SW-OT17B**

- A fiber optic otoscope(OT10)
- A direct ophthalmoscope(OP10)
- 2.4mm/3mm/4mm/5mm tips tow groups
- Soft nylon bag
- 一个光纤检耳镜 (OT10)
- 一个直接检眼镜 (OP10)
- 两组2.2mm/3mm/4mm/5mm 耳套
- 尼龙布袋

Package content

Mode	Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
SW-OT17SA	193x102x32	400x215x195	20	5.2/4.8
SW-OT17	190x110x40	400x240x230	20	7.5/7.0
SW-OT17C	190x110x48	400x215x265	20	6.5/6.0
SW-OT17B	190x110x40	400x240x230	20	8.4/7.9

Gross and net weight(Kgs/Ctn)Battery not included

**SW-OT16LK**

- A deluxe LCD screen handle
- 2.4mm/3mm/4mm/5mm tips of tow groups
- A head of fiber optic otoscope
- A head of laryngeal mirror
- 一个豪华液晶屏手柄
- 两组2.4mm/3mm/4mm/5mm 耳套
- 一个光纤检耳镜头部
- 一个口腔镜头部

**SW-OT15LK**

- A deluxe LCD screen handle
- A head of direct ophthalmoscope
- A head of laryngeal mirror
- 一个豪华液晶屏手柄
- 一个直接检眼镜头部
- 一个口腔镜头部

**SW-OT17L**

- A deluxe LCD screen handle
- 2.4mm/3mm/4mm/5mm tips of tow groups
- A head of fiber optic otoscope
- A head of direct ophthalmoscope
- 一个豪华液晶屏手柄
- 两组2.4mm/3mm/4mm/5mm 耳套
- 一个光纤检耳镜头部
- 一个直接检眼镜头部

**SW-OT17LK**

- A deluxe LCD screen handle
- 2.4mm/3mm/4mm/5mm tips of tow groups
- A head of fiber optic otoscope
- A head of direct ophthalmoscope
- A head of laryngeal mirror
- 一个豪华液晶屏手柄
- 两组2.4mm/3mm/4mm/5mm 耳套
- 一光纤检耳镜头部
- 一个直接检眼镜头部
- 一个口腔镜头部



- The handle of the above series of products can be selected as a standard type of handle
- 上面套装系列产品的手柄可选配为标准型手柄

Package content

Mode	Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
SW-OT16LK	190x110x40	400x240x230	20	7.6/7.1
SW-OT15LK	190x110x40	400x240x230	20	8.1/7.6
SW-OT17L	190x110x40	400x240x230	20	7.7/7.2
SW-OT17LK	190x110x40	400x240x230	20	8.9/8.4

Gross and net weight(Kgs/Ctn)Battery not included

**SW-OT18LKY**

- A deluxe LCD screen handle
- 2.4mm/3mm/4mm/5mm tips of tow groups
- A head of fiber optic otoscope
- A head of laryngeal mirror and reflector lens
- A head of tongue depressor holde
- 5 pcs tongue depressor
- 一个豪华液晶屏手柄
- 两组2.4mm/3mm/4mm/5mm 耳套
- 一个光纤检耳镜头部
- 一个口腔镜头部和反光镜头
- 一个压舌板器头部
- 5片压舌板

**SW-OT18LKE**

- A deluxe LCD screen handle
- 2.2mm/3mm/4mm/5mm/9mm tips is a groups
- A head of direct ophthalmoscope
- A head of fiber optic otoscope
- A head of nasal speculum
- 一个豪华液晶屏手柄
- 一组2.2mm/3mm/4mm/5mm /9mm耳套
- 一个直接检眼镜头部
- 一个口腔镜头部和反光镜头
- 一个鼻腔镜头部

**SW-OT18YKE**

- A deluxe LCD screen handle
- 2.2mm/3mm/4mm/5mm/9mm tips is a groups
- A head of laryngeal mirror and reflector lens
- A head of nasal speculum
- A head of tongue depressor holder
- 5 pcs tongue depressor
- 一个豪华液晶屏手柄
- 一组2.2mm/3mm/4mm/5mm /9mm耳套
- 一个口腔镜头部和反光镜头
- 一个鼻腔镜头部
- 一个压舌板器头部
- 5片压舌板

**SW-OT15LKEY**

- A deluxe LCD screen handle
- 2.2mm/3mm/4mm/5mm/9mm tips is a groups
- A head of direct ophthalmoscope
- A head of laryngeal mirror and reflector lens
- A head of nasal speculum
- A head of tongue depressor holder
- 5 pcs tongue depressor
- 一个豪华液晶屏手柄
- 一组2.2mm/3mm/4mm/5mm /9mm耳套
- 一个直接检眼镜头部
- 一个口腔镜头部和反光镜头
- 一个鼻腔镜头部
- 一个压舌板器头部
- 5片压舌板

Package content

Mode	Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
SW-OT18LKY	275x160x56	340x295x310	10	8.4/7.7
SW-OT18LKE	275x160x56	340x295x310	10	8.4/7.7
SW-OT18YKE	275x160x56	340x295x310	10	9.0/8.3
SW-OT18LKEY	275x160x56	340x295x310	10	9.0/8.3

Gross and net weight(Kgs/Ctn)Battery not included



- The handle of the above series of products can be selected as a standard type of handle
- 上面套装系列产品的手柄可选配为标准型手柄

**SW-OT16C**

- A rechargeable fiber optic otoscope
- 一个可充电光纤检耳镜
- 2.2mm/3mm/4mm/5mm/9mm tips is a groups
- 一组2.2mm/3mm/4mm/5mm/9mm耳套
- AA Ni-Mh2.4V rechargeable battery
- 2.4V AA镍氢充电电池
- Single-charger
- 单孔充电座

- Search bar and lens are removable
- 可拆卸的探索杆和镜头
- Straight rod or Curved rod are optinal
- 弯杆直杆可选择

**SW-OT15C**

- 5 different aperturus
- 五种光阑
- 19 kinds of adjustable corrective lenses,
- Range of lenses : -20d ~ +20d
- 19种屈光镜片, 范围: -20D ~ +20D

- 2x optical glass observation window lens
- 2x光学玻璃观察窗镜片
- It can been upgrade to the HS-OT10H
- 可升级为HS-OT10H手术耳镜
- 5 pcs 9mm nose tips
- 5个9mm窥鼻器

**SW-OT18Y3C**

- Patented tongue-depressor ejectable device
- 拥有专利的压舌板弹出装置

- 3200k and 6000k color temperature LED bulb and XHL Xenon Halogen bulb to optiona
- 3200k和6000k色温LED灯和XHL卤素灯可选

**SW-OT18K3C****SW-OT18E3C****Package content**

Mode	Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
SW-OT16C	206x173x88	460x390x330	12	11.6/10.6
SW-OT18K3C	206x173x88	460x390x330	12	11.5/10.5
SW-OT15C	206x173x88	460x390x330	12	11.0/10.0
SW-OT18E3C	206x173x88	460x390x330	12	11.8/10.8
SW-OT18Y3C	206x173x88	460x390x330	12	11.9/10.9

Gross and net weight(Kgs/Ctn)Battery not included



HS-OT10L-C

- A rechargeable fiber optic otoscope
一个可充电光纤检耳镜
- 2.2mm/3mm/4mm/5mm/9mm tips is a groups
一组2.2mm/3mm/4mm/5mm/9mm耳套
- AA Ni-Mh2.4V rechargeable battery
2.4V Aa镍氢充电电池
- Single-charger
单孔充电座



HS-OP10L-C

- 5 different aperturus
五种光阑
- 19 kinds of adjustable corrective lenses,
Range of lenses : -20d~+20d
19种屈光镜片, 范围:-20D~+20D



HS-OT10Y-C

- Patented tongue-depressor ejectable device
拥有专利的压舌板弹出装置



Start charging



Finish charging

Package content

Mode	Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
SW-OT16LC	206x173x88	460x390x330	12	11.6/10.6
SW-OT18KC	206x173x88	460x390x330	12	11.5/10.5
SW-OT15LC	206x173x88	460x390x330	12	11.0/10.0
SW-OT18EC	206x173x88	460x390x330	12	11.8/10.8
SW-OT18YC	206x173x88	460x390x330	12	11.9/10.9

Gross and net weight(Kgs/Ctn)Battery not included



HS-OT10K-C

- Search bar and lens are removable
可拆卸的探索杆和镜头
- Straight rod or Curved rod are optional
弯杆直杆可选择



HS-OT10E-C



SW-OT17C2

- A rechargeable fiber optic otoscope
一个可充电光纤检耳镜
- A rechargeable ophthalmoscope
一个可充电的检眼镜
- 2.4mm/3mm/4mm/5mm tips of two groups
两组2.4mm/3mm/4mm/5mm耳套
- One insufflation bulb
一个鼓气球
- Double charger
双孔充电座
- Gift box of set package
礼盒套装



SW-OT18H3C2

- A rechargeable fiber optic otoscope
一个可充电光纤检耳镜
- A operating otoscope of chargeable
一个可充电手术耳镜
- 2.2mm/3mm/4mm/5mm/9mm tips is a groups
一组2.2mm/3mm/4mm/5mm/9mm耳套



SW-OT18Y3C2

- A rechargeable fiber optic otoscope
一个可充电光纤检耳镜
- A tongue depressor of chargeable
一个可充电压舌板

- 3200k and 6000k color temperature LED bulb and XHL Xenon Halogen bulb to optiona
3200k和6000k色温LED灯和XHL卤素灯可选



SW-OT183C2

- A rechargeable fiber optic otoscope
一个可充电光纤检耳镜
- A laryngeal mirror of chargeable
一个可充电口腔镜



SW-OT18E3C2

- A rechargeable fiber optic otoscope
一个可充电光纤检耳镜
- A nasal speculum of chargeable
一个可充电鼻腔镜
- 5 pcs 9mm nose tips
5个9mm窥鼻器

Package content

Mode	Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
SW-OT17C2	206x173x88	460x390x330	12	14.6/13.6
SW-OT183C2	206x173x88	460x390x330	12	14.3/13.3
SW-OT18H3C2	206x173x88	460x390x330	12	14.7/13.7
SW-OT18E3C2	206x173x88	460x390x330	12	15.0/14.0
SW-OT18Y3C2	206x173x88	460x390x330	12	14.5/13.5

Gross and net weight(Kgs/Ctn)Battery not included



SW-OT17LCS

- A rechargeable fiber optic otoscope
一个可充电光纤检耳镜
- A rechargeable ophthalmoscope
一个可充电的检眼镜
- 2.4mm/3mm/4mm/5mm tips of two groups
两组2.4mm/3mm/4mm/5mm耳套
- One insufflation bulb
一个鼓气球
- Double charger
双孔充电座
- Gift box of set package
礼盒套装



SW-OT18HC2

- A rechargeable fiber optic otoscope
一个可充电光纤检耳镜
- A laryngeal mirror of chargeable
一个可充电口腔镜

- A rechargeable fiber optic otoscope
一个可充电光纤检耳镜
- A operating otoscope of chargeable
一个可充电手术耳镜
- 2.2mm/3mm/4mm/5mm/9mm tips is a groups
一组2.2mm/3mm/4mm/5mm/9mm耳套



SW-OT18YC2

- A rechargeable fiber optic otoscope
一个可充电光纤检耳镜
- A nasal speculum of chargeable
一个可充电鼻腔镜
- 5 pcs 9mm nose tips
5个9mm窥鼻器

- A rechargeable fiber optic otoscope
一个可充电光纤检耳镜
- A tongue depressor of chargeable
一个可充电压舌板

Start charging Finish charging



Patent protection

Package content

Mode	Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
SW-OT17LCS	206x173x88	460x390x330	12	14.6/13.6
SW-OT18KC2	206x173x88	460x390x330	12	14.3/13.3
SW-OT18HC2	206x173x88	460x390x330	12	14.7/13.7
SW-OT18EC2	206x173x88	460x390x330	12	15.0/14.0
SW-OT18YC2	206x173x88	460x390x330	12	14.5/13.5

Gross and net weight(Kgs/Ctn)Battery not included



SW-OT18KC2



SW-OT18EC2

SW-ENT100 Integrated Diagnostic System



MAIN POWER
• INPUT: AC220V/110V
• OUTPUT: DC3.5V



[M1] FIBER OPTIC OTOSCOPE

- Lightweight design with high performance
- Tunable Light Source
- 3X Magnification
- Fiber optic illumination
- Xenon Halogen Technology/LED Technology (optional)
- Remark:with insufflation bulb

[M1] 光纤耳镜

- 轻量化设计, 高性能
- 光源柔和
- 3X大视窗
- 光纤照明
- 卤素/LED灯 (可选)
- 气囊球



[M3] INFRARED THERMOMETER

- Forehead thermometer, 1 second measurement
- Multiple use (Wide Range Measurement)
- °C<->°F Switchable
- 30 Data memory
- Last read memory
- Fever alarm

[M3] 额温枪

- 1秒测量
- 多样使用 (测试范围广)
- °C-°F可切换
- 30天内存
- 记忆最后读数
- 发热报警装置

Insufflation Bulb
气囊球

- M1&M2 handle can be equipped with other models: W1、W2、W3、W4、W5 heads
- M1和M2手柄可配其他型号: W1、W2、W3、W4、W5的头



[M2] DIRECT OPHTHALMOSCOPE

- Lightweight design and easy to operation
- Tunable Light Source
- 6 apertures with red-free and blue filter
- 29 kinds of adjustable corrective Range of lenses :-35D~+40D
- Xenon Halogen Technology

[M2] 直接检眼镜

- 轻量化设计, 简易操作
- 光源柔和
- 6种光阑和无赤滤光、钴蓝滤光
- 29种屈光透镜, 屈光范围: -35D~+40D
- 卤素/LED灯 (可选)



[M4] SPHYGMOMANOMETER

- Ergonomics design
- The Square shaped dial 145X145mm
- Back box is easy to hold the cuff
- High quality Sphygmomanometer

[M4] 血压表

- 人体工学设计
- 145mm*145mm的正方形设计
- 后背盒子可方便的放置产品袖口绑带
- 高质量血压表

[M1.1] TIP DISPENSER
(WITH DISPOSABLE TIPS)

- Tip Size: 2.75mm dia. (2 Sets)
- 4.25mm dia. (2 Sets)
- Each Set: 15 pieces

[M1.1] 耳套 (一次性耳套)

- 耳套尺寸: 2.75mm (2组)
- 4.25mm (2组)
- 每种15个



W1



W2



W3




W4



W5

SW-100ENTW



[W1] Fiber Optic Otoscope

- 3X magnification for the large size viewing window
- Fiber optic illumination
- The socket of ear insufflation bulb
- Supplied by 2xAAA alkaline batteries (not included)
- The LCD screen handle can adjust the brightness of the light source
- Optical SMT LED technology, CRI≥90%
- Constant current IC circuit, maintain the optimal lighting effect
- Modular stainless steel bracket is convenient and durable



- 3倍放大的观察视窗
- 光纤照明
- 气囊球接口
- 2XAAA碱性电池 (不包含)
- 液晶屏手柄, 可以调节光源亮度
- Optical SMT LED 技术, 显指≥90%
- 内置恒流稳压芯片, 光效稳定
- 模块化不锈钢支架, 方便耐用



[W1.1] Tip Dispenser

- Tip size: 2.4mm/3mm/4mm/5mm dia
- Each model: 12 pieces
- Single tube arrangement easy to extract
- module bracket save space

- 耳套尺寸2.4mm/3mm/4mm/5mm
- 每个型号12个
- 模块支架, 节省空间
- 单管排列, 抽取方便



[W4] Nasal Speculum



[W2] Direct Ophthalmoscope

- 5 different apertures
- 19 kinds of adjustable corrective lenses, Range of lenses: -20D~+20D
- The LCD screen handle can adjust the brightness of the light source
- CRI≥90%, it is more color reducibility
- Optical SMT LED technology, low power consumption and long life

- 五种光阑
- 19种屈光透镜, 屈光度范围: -20D~+20D
- 液晶屏手柄, 可以调节光源亮度
- 显指≥90% 色彩还原性强
- Optical SMT LED技术, 低功耗、长寿命



[W3] Laryngeal Mirror

- Search rod and lens are removable to convenient detoxification
- Anti-fog lens the clear observation
- Straight rod or Curved rod are optional
- Optional disposable laryngeal mirror head for easy installation and disassembly


- 可拆卸的探索杆及镜头, 方便消毒
- 防雾镜片, 观察清晰
- 直杆, 弯杆可选
- 可选配一次性镜头, 便于安装和拆卸



[W5] Tongue Depressor Holder

- It is more convenient to use with tongue-depressor ejectable device
- Plastic and standard wooden spatula can be employed universally

- 压舌板弹出装置, 使用更安全
- 塑料压舌板和一次木质压舌板可通用



[W8] Sphygmomanometer

- Blood pressure classification indication
- Irregular heart-beat detector
- 3 sets of alarms
- Ambient temperature display during sleep mode
- Date & Time display

- 血压分级指示
- 心率不全检测
- 3组报警提示
- 休眠模式下显示环境温度
- 日期和时间显示



[W7] Infrared Thermometer

- Fahrenheit/centigrade switchable
- 2 seconds fast reading
- Accuracy degree: ±0.1°C
- Heat alarm function
- 2x1.5 AAA by alkaline battery

- 华氏/摄氏度可切换
- 2秒快速读取
- 精确度: ±0.1°C
- 发热报警
- 使用2*1.5V AAA碱性电池



[W10] Blood Glucose Monitoring System

- 5-sec display testing result
- 0.5ul blood required tiny quantity
- Storable 900 test results

- 5秒获取测试结果
- 0.5ul微量血液
- 存储900条测试结果



[W9] Stethoscope

- Silver dual-membrane chestpiece
- Excellent dual-membrane chestpiece, little acoustic attenuation
- Standard dual PVC tube, good resistance to ambient interference

- 银质双膜片头
- 良好的导声性能和小的声音衰减
- 标准的PVC双管, 抗环境干扰强



ENT Diagnostic Series

W3



- Laryngeal mirror
- 口腔镜头部

W4



- Nasal speculum
- 鼻腔镜头部

W5



- Tongue-depressor holder
- 压舌板器头部



ENT


SW-OT18M



Lr4s LED technical parameters

1. Power: 0.2w
2. Colour temperature: 3000K / 5000K
3. Luminance: 4.5Lm
4. Illuminance: 21000Lux, light intensity 0.25m/5cd, Stereoscopic 30 degrees uniform
5. Battery: 1pc disposable battery AAA/1.5V, endurance 4hours.
1pc rechargeable battery AAA/1.2V

LED luminance of high quality curve graph




1. Otoloscope



B. Adult/Child size earpick with light

2. Conchoscope. Metal: diameter is 5.0mm
Light bending rod

3. Mouth speculum



4. With tongue depressor




5. Eye examination

- A. The flashlight is directed at 30 degrees and the low light is 0.25m/5cd, and the surface of the eye membrane can be clearly inspected.



- B. The flashlight with a soft light rod (magnetic), and the light rod is against the inner layer of the eyelid to remove foreign bodies, debris and small particles.

Optical Fiber
oral cavity mirrorTongue
depressorNasal
optical fiberMagnetic
light stickEar spoon
light stick

Otoscope

TOURNIQUET

SW-GT26



Packing: 1pc/bag, 100pcs/box,
500pcs/CTN,
50x40x27cm 14kgs

SW-GT03



Packing: 1pc/bag, 25pcs/box,
500pcsc/CTN,
50x40x27cm, 14kgs

SW-GT01



Packing: 1pc/bag, 25pcs/box,
500pcs/ctn,
50x40x27cm, 19kgs

SW-GT06



Packing: 1pc/bag, 25pcs/box,
500pcs/ctn,
50x40x27cm, 19kgs

SW-GT18



Packing: 1pc/bag, 25pcs/box,
500pcs/ctn,
50x40x27cm, 16kgs

SW-GT19



Packing: 1pc/bag, 25pcs/box,
500pcs/ctn,
50x40x27cm, 17kgs

SW-GT14



Packing: 1pc/bag, 25pcs/box,
500pcs/ctn,
50x40x27cm, 19kgs

SW-GT04



Packing: 1pc/bag, 50pcs/box,
1000pcs/ctn,
50x40x27cm, 15kgs

SW-GT08



Disposable non-latex tourniquet
Material: TPE size: 1"x18"
Packing: tiled, 250pcs/box, 1000pcs/ctn,
50x19x16cm, 6.8kgs
roll, 500pcs/bag, 1000pcs/ctn,
37x24.5x18.5cm, 6.4kgs

PENLIGHT

SW-PL07



SW-OT01

1. with 5pcs button batteries, LED bulb
 2. with 2pcs AAA batteries, Halogen bulb
- Packing: 1pc/box, 100pcs/ inner box, 500pcs/ctn,
78x20x30cm 11kgs 15kgs
with tongue depressor holder

SW-PL01



- with 5pcs button batteries, LED bulb
packing: 6pcs/card bag, 60pcs/box, 600pcs/ctn,
48x18x27cm 10.5kgs
with 2pcs AAA batteries , Halogen bulb, 15kgs

SW-PL78



- Stainless steel shell LED bulb
Super operation time:10000 hours 2pcs AAA batteries

SW-OT03



- Mini otoscope
1. with 5pcs button batteries, LED bulb
2. with 2pcs AAA batteries, Halogen bulb
Packing: 1pc/box, 100pcs/ inner box, 500pcs/ctn,
78x20x30cm 11kgs 15kgs

PENLIGHT

SW-G81



Depressor holder Available mode(SW-PL07)

SW-PL45



1. With LED bulb
 2. With Halogen bulb
- Packing: 1pc/box, 100pcs/inner box, 33x34x48cm
500pcs 10kgs

SW-PL48



- LED light
Packing: 1pc/box, 100pcs/inner box,
33x34x48cm 500pcs 10kgs

SW-PL49



- With Halogen bulb
Packing: 1pc/box, 100pcs/inner box,
33x34x48cm 500pcs 10kgs

NEUROLOGICAL HAMMER

SW-H01



Packing: 1pcs/box 100pcs/CTN, 39x33x21.5cm, 11kgs

Zinc alloy handle; PVC head; Diagnostic Pin can be screwed off conveniently in the head, while diagnostic brush can be taken out in the end.

SW-H02



Packing: 1pcs/box 100pcs/CTN, 48.5x45x27cm, 10kgs

Zinc alloy handle; PVC head; Color is optional.

SW-H03



Packing: 1pcs/box 100pcs/CTN, 40x19x23cm, 5kgs

Zinc alloy handle; PVC head; Copper head.

SW-H04



Packing: 1pcs/box 100pcs/CTN, 40x29.5x23cm, 7kgs

Zinc alloy handle; PVC head; Color is optional.

NEUROLOGICAL HAMMER

SW-H05



Packing: 1pcs/box 100pcs/CTN, 49.5x48x17cm, 11.6kgs

Copper handle with rubber painting; PVC head; Color is optional.

SW-H06



Packing: 1pcs/box 100pcs/CTN, 33x23x23cm, 9kgs

Stainless steel handle; PVC head; Color is optional.

SW-H07



Packing: 1pcs/blister, 100pcs/CTN, 55.5x34x17cm, 11kgs

Zinc alloy handle; Zinc alloy+silicone head; Color is optional.

SW-H10



Packing: 1pcs/box 150pcs/CTN, 40x39x20cm, 14kgs

NEUROLOGICAL HAMMER

SW-H16



Packing: 1pcs/box 50pcs/CTN, 47.5x25x31cm

Zinc alloy; Silicone head; Color is optional.

SW-H41



Packing: 1pcs/box, 100pcs/CTN, 44x27x32cm, 10kgs

Zinc alloy+ABS handle; Silicone head.

SW-H42



Packing: 1pcs/box, 100pcs/CTN, 35x27.5x25cm, 5kgs

ABS handle; silicone head; Color is optional.

SW-H45



Packing: 1pcs/blister 100pcs/CTN, 55.5x34x17cm, 7.3kgs

ABS handle; PVC head; Color is optional.

NEUROLOGICAL HAMMER

SW-H47



Packing: 1pcs/blister 50pcs/CTN, 42x41x21cm, 7.8kgs

ABS handle; Zinc alloy + silicone head; Color is optional.

SW-H46



Packing: 1pcs/box 100pcs/CTN, 40.5x30.5x26cm, 9kgs

ABS handle with rubber painting. PVC head; Color is optional.

SW-H48



Packing: 1pcs/box, 100pcs/CTN, 50x24x27cm, 7kgs

ABS handle; Silicone head; Color is optional.

SW-H49



Packing: 1pcs/blister, 50pcs/CTN, 42x41x21cm, 9kgs

Nylon handle; Zinc alloy + silicone head; Color is optional.

NEUROLOGICAL HAMMER

SW-H50



Packing: 1pcs/blister 50pcs/CTN, 42x41x21cm, 7kgs

Zinc alloy handle; Silicone head; Color is optional.

SW-H58



Packing: 1pcs/blister 100pcs/CTN, 37x31x34cm, 10kgs

ABS handle; Silicone head; Color is optional.

SW-H59



Packing: 1pcs/blister 50pcs/CTN, 42x41x21cm, 8.5kgs

ABS handle; Zinc alloy + silicone head; Color is optional.

SW-H60



Packing: 1pcs/blister, 50pcs/CTN, 42x41x21cm, 6.5kgs

Nylon handle; Zinc alloy + silicone head; Color is optional.

MIX SET

SW-GP06	Detailed				
					
	Diagnostic kit Penlight + Otoscope head + Tips + Spatula + 2xAAA batteries 套装包含 手电筒 + 检耳镜头 + 5个耳套 + 压舌板 + 2节7号电池				
Packing Bill 包装	Box size 中盒尺寸	Carton size 外箱尺寸	Pcs/Ctn 数量	Kgs/Ctn 重量	
Plastic box 塑料盒	210x110x30 (mm)	570x440x310 (mm)	100	15.6	

SW-OT11	Detailed				
					
	Diagnostic kit Penlight + Otoscope head + Tips + 2xAAA batteries. 套装包含 手电筒 + 检耳镜头 + 5个耳套 + 2节7号电池				
Packing Bill 包装	Box size 中盒尺寸	Carton size 外箱尺寸	Pcs/Ctn 数量	Kgs/Ctn 重量	
Painted Inner box	200x85x30 (mm)	500x400x270 (mm)	100	11.5	

SW-GP09B	Detailed				
					
	Diagnostic kit Penlight + Otoscope head + Spatula clip + Spatula + Blue filter + Ear pick + 2xAAA batteries 套装包含 手电筒 + 检耳镜头 + 压舌板夹 + 压舌板 + 蓝色过滤器 + 耳勺 + 2节7号电池				
Packing Bill 包装	Box size 中盒尺寸	Carton size 外箱尺寸	Pcs/Ctn 数量	Kgs/Ctn 重量	
Plastic box 塑料盒	188x100x25 (mm)	515x390x265 (mm)	100	14	


SW-GP09C	Detailed				
					
	Diagnostic kit Penlight + Otoscope head + Tips + Blue filter + Spatula clip + Spatula + 2xAAA batteries 套装包含 手电筒 + 检耳镜头 + 3个耳套 + 蓝色过滤器 + 压舌板夹 + 压舌板 + 2节7号电池				
Packing Bill 包装	Box size 中盒尺寸	Carton size 外箱尺寸	Pcs/Ctn 数量	Kgs/Ctn 重量	
Plastic box 塑料盒	188x100x25 (mm)	515x390x265 (mm)	100	14	

MIX SET

SW-GP20			Detailed		
			Diagnostic kit Reflex hammer + Otoscope head + Penlight + 2xAAA batteries		
	Packing Bill 包装	Box size 中盒尺寸	Carton size 外箱尺寸	Pcs/Ctn 数量	Kgs/Ctn 重量
	Plastic box 塑料盒	215x120x30 (mm)	460x260x305 (mm)	40	13


SW-GP09D			Detailed		
			Diagnostic kit Penlight + Otoscope head + Tips + 2xAA batteries. 套装包含 手电筒 + 检耳镜头 + 3个耳套 + 2节5号电池		
	Packing Bill 包装	Box size 中盒尺寸	Carton size 外箱尺寸	Pcs/Ctn 数量	Kgs/Ctn 重量
	Plastic box 塑料盒	190x70x30 (mm)	400x350x330 (mm)	100	16.5

SW-GP01			Detailed		
			Diagnostic kit 5pcs hammers mix set		
	Packing Bill 包装	Box size 中盒尺寸	Carton size 外箱尺寸	Pcs/Ctn 数量	Kgs/Ctn 重量
	Painted Inner Box	345x130x30 (mm)	570x365x300 (mm)	40	23

SW-GP02			Detailed		
			Diagnostic kit Pinwheel hammer + Reflex hammer + Penlight		
	Packing Bill 包装	Box size 中盒尺寸	Carton size 外箱尺寸	Pcs/Ctn 数量	Kgs/Ctn 重量
	Painted Inner Box	345x130x30 (mm)	570x365x300 (mm)	40	16

MONOFILAMENT

SW-G307



Packing: 1pc/bag, 500pcs/CTN,
80x24x23cm 10kgs
Five probes; ABS shell;
Color is optional.

SW-G30



Packing: 1pc/box, 600pcs/CTN,
57x38x18.5cm 12.5kgs
PP material.

SW-G301



Packing: 1pc/bag, 500pcs/CTN,
55x33.5x30cm 10kgs
Five probes; ABS shell;
Color is optional.

SW-G302



Packing: 1pc/bag, 500pcs/CTN,
42x30x40cm 9kgs
ABS material.

SW-G303



Packing: 1pc/bag, 500pcs/CTN,
37x23x43cm 9kgs
ABS material.

SW-G304



Packing: 1pc/bag, 2000pcs/CTN,
48x30x33cm 19.5kgs
ABS material.

SW-G306



Packing: 1pc/bag, 2000pcs/CTN,
48x30x33cm 19.5kgs
ABS material.

OXIMETER



SW-FO1A3

Features

- ▲ Can measure SpO2 and Pulse Rate accurately
- ▲ SpO2 and Pulse Rate display, Pulse Rate waveform and bar graph display
- ▲ Battery voltage low indication
- ▲ OLED display(two color), Four directions and six modes
- ▲ Screen brightness adjustable
- ▲ Shut off automatically when no signal
- ▲ 6 Color available: Grey, Green, Purple, Blue, Pink, Brown

SpO2	Measurement range: 70-99% Resolution: ±1% Accuracy: ±2% (70%-99%), unspecified (<70%)
Pulse rate	Measurement range: 30-240 bpm Resolution: ±1% Accuracy: ±2bpm or ±2% (select target) Low Perfusion: ≤0.4%
Power	1.5V (AAA size) alkaline battery x 2 Supply voltage: 2.6-3.6V
Working current	≤30mA
Automatic power-off	Automatically power off when no signal in the oximeter for more than 8 seconds
Dimension & Weight	62 (L) x 90 (W) x 40 (H) mm 65g (without batteries)

SW-FO1B1

Features

- ▲ Can measure SpO2 and Pulse Rate accurately
- ▲ SpO2 and Pulse Rate display, Pulse Rate waveform and bar graph display
- ▲ Battery voltage low indication
- ▲ OLED display(two color), Four directions and six modes
- ▲ Screen brightness adjustable
- ▲ Shut off automatically when no signal
- ▲ 5 Color available: Green, Red, Yellow, Black, Blue, White



SpO2	Measurement range: 70-99% Resolution: ±1% Accuracy: ±2% (70%-99%), unspecified (<70%)
Pulse rate	Measurement range: 30-240 bpm Resolution: ±1% Accuracy: ±2bpm or ±2% (select target) Low Perfusion: ≤0.4%
Power	1.5V (AAA size) alkaline battery x 2 Supply voltage: 2.6-3.6V
Working current	≤30mA
Automatic power-off	Automatically power off when no signal in the oximeter for more than 8 seconds
Dimension & Weight	62 (L) x 90 (W) x 40 (H) mm 65g (without batteries)

SW-FO1B2

Features

- ▲ Can measure SpO2 and Pulse Rate accurately
- ▲ SpO2 and Pulse Rate display, Pulse Rate waveform and bar graph display
- ▲ Battery voltage low indication
- ▲ OLED display(two color), Four directions and six modes
- ▲ Screen brightness adjustable
- ▲ Shut off automatically when no signal
- ▲ 5 Color available: Green, Red, Yellow, Blue, White



SpO2	Measurement range: 70-99% Resolution: ±1% Accuracy: ±2% (70%-99%), unspecified (<70%)
Pulse rate	Measurement range: 30-240 bpm Resolution: ±1% Accuracy: ±2bpm or ±2% (select target) Low Perfusion: ≤0.4%
Power	1.5V (AAA size) alkaline battery x 2 Supply voltage: 2.6-3.6V
Working current	≤30mA
Automatic power-off	Automatically power off when no signal in the oximeter for more than 8 seconds
Dimension & Weight	60 (L) x 90 (W) x 38 (H) mm 54g (without batteries)

SW-FO1C1

Features

- ▲ Can measure SpO2 and Pulse Rate accurately
- ▲ SpO2 and Pulse Rate display, Pulse Rate waveform and bar graph display
- ▲ Battery voltage low indication
- ▲ OLED display(two color), Four directions and six modes
- ▲ Screen brightness adjustable
- ▲ Shut off automatically when no signal
- ▲ 4Color available: Red, Orange, Blue, White



SpO2	Measurement range: 70-99% Resolution: ±1% Accuracy: ±2% (70%-99%), unspecified (<70%)
Pulse rate	Measurement range: 30-240 bpm Resolution: ±1% Accuracy: ±2bpm or ±2% (select target) Low Perfusion: ≤0.4%
Power	1.5V (AAA size) alkaline battery x 2 Supply voltage: 2.6-3.6V
Working current	≤30mA
Automatic power-off	Automatically power off when no signal in the oximeter for more than 8 seconds
Dimension & Weight	60 (L) x 92 (W) x 42 (H) mm 65g (without batteries)

OXIMETER



SW-FO1H1

Features

- ▲ Can measure SpO2 and Pulse Rate accurately
- ▲ SpO2 and Pulse Rate display, Pulse Rate waveform and bar graph display
- ▲ Battery voltage low indication
- ▲ OLED display(two color), Four directions and six modes
- ▲ Screen brightness adjustable
- ▲ Shut off automatically when no signal
- ▲ 6 Color available: Gray, Apricot, Blue, Green, Pink, Purple

SpO2	Measurement range: 70-99% Resolution: ±1% Accuracy: ±2% (70%-99%), unspecified (<70%)
Pulse rate	Measurement range: 30-240 bpm Resolution: ±1% Accuracy: ±2bpm or ±2% (select target) Low Perfusion: ≤0.4%
Power	1.5V (AAA size) alkaline battery x 2 Supply voltage: 2.6-3.6V
Working current	≤30mA
Automatic power-off	Automatically power off when no signal in the oximeter for more than 8 seconds
Dimension & Weight	100 (L) x 64 (W) x 37 (H) mm 65g (without batteries)

SW-FO1J1

Features

- ▲ Can measure SpO2 and Pulse Rate accurately
- ▲ SpO2 and Pulse Rate display, Pulse Rate waveform and bar graph display
- ▲ Battery voltage low indication
- ▲ OLED display(two color), Four directions and six modes
- ▲ Screen brightness adjustable
- ▲ Shut off automatically when no signal
- ▲ 2Color available: Blue and Yellow



SpO2	Measurement range: 70-99% Resolution: ±1% Accuracy: ±2% (70%-99%), unspecified (<70%)
Pulse rate	Measurement range: 30-240 bpm Resolution: ±1% Accuracy: ±2bpm or ±2% (select target) Low Perfusion: ≤0.4%
Power	1.5V (AAA size) alkaline battery x 2 Supply voltage: 2.6-3.6V
Working current	≤30mA
Automatic power-off	Automatically power off when no signal in the oximeter for more than 8 seconds
Dimension & Weight	80 (L) x 60 (W) x 52 (H) mm 70g (without batteries)

SW-FO1JO

Features

- ▲ Can measure SpO2 and Pulse Rate accurately
- ▲ Battery voltage low indication
- ▲ Easy to read LED display
- ▲ Shut off automatically when no signal
- ▲ 2 Color available: Blue and Yellow




SpO2	Measurement range: 70-99% Resolution: ±1% Accuracy: ±2% (70%-99%), unspecified (<70%)
Pulse rate	Measurement range: 30-240 bpm Resolution: ±1% Accuracy: ±2bpm or ±2% (select target) Low Perfusion: ≤0.4%
Power	1.5V (AAA size) alkaline battery x 2 Supply voltage: 2.6-3.6V
Working current	≤30mA
Automatic power-off	Automatically power off when no signal in the oximeter for more than 8 seconds
Dimension & Weight	80 (L) x 62 (W) x 42 (H) mm 65g (without batteries)

SW-FO80A

Features

- ▲ Can measure SpO2 and Pulse Rate accurately
- ▲ SpO2 and Pulse Rate display, Pulse Rate waveform and bar graph display
- ▲ Battery voltage low indication
- ▲ OLED display(two color), Four directions and six modes
- ▲ Screen brightness adjustable
- ▲ Shut off automatically when no signal
- ▲ 5 Color available: Grey, Green, Purple, Blue, Pink



OXIMETER

SW-FO82

Fingertip Pulse Oximeter

Display: LED

Parameters: SpO2, PR, Waveform, Pulse bar

Features:

- ▲ Gravity function, Automatic rotating read
- ▲ PI-Perfusion Index indication
- ▲ Dual color OLED display SpO2, PR, Waveform, Pulse bar
- ▲ 4-direction & 6 mode display provide convenient reading
- ▲ Setting alarm range of SpO2 and pulse rate
- ▲ Menu function setting (Beep sounds, etc)
- ▲ 2pcs AAA alkaline batteries, automatically power off



SW-FO84

OLED fingertip oximeter

- ▲ Dual color OLED display SpO2, PR, waveform, pulse bar
- ▲ 4 direction, 6 mode display provide convenient reading
- ▲ PI-perfusion index indication
- ▲ Setting alarm range of Spo2 and pulse rate
- ▲ Menu function setting (beep sound, etc)
- ▲ 2 x AAA alkaline batteries
- ▲ Gravity function, automatic rotating read (Option)
- ▲ 57.5x35.8x30mm, 27g



SW-FOC1

Fingertips Oximeter

- ▲ LED full screen display
- ▲ 2 x AAA 1.5V alkaline batteries
- ▲ 8seconds automatic shutdown
- ▲ Color: black, white
- ▲ Size: 58x36x33mm



OXIMETER

SW-FOC1B

OLED Fingertips Oximeter\

- ▲ Dual color OLED full screen display
- ▲ 4 display directions , 6 display styles
- ▲ 2 x AAA 1.5V alkaline batteries
- ▲ Brightness can adjustable from 1 to 4
- ▲ 8 seconds automatic shutdown
- ▲ Color: black, white
- ▲ Size 58x36x33cm



SW-FOT1

Fingertip Oximeter

- ▲ Dual color OLED full screen display
- ▲ 4 display directions, 6 display styles
- ▲ 2 x AAA alkaline batteries
- ▲ Brightness can adjustable from 1 to 4
- ▲ 8 seconds automatic shutdown
- ▲ Size: 59x33x30mm



SW-FOT2

Fingertip Oximeter

- ▲ Dual color OLED full screen display
- ▲ 4 display directions, 6 display styles
- ▲ 2 x AAA alkaline batteries
- ▲ Brightness can adjustable from 1 to 4
- ▲ 8 seconds automatic shutdown
- ▲ Size: 59x33x30mm



MEDICAL GIFTS

SW-G06

High measurement



SW-G37
Pen holder

SW-G06B

High measurement

Item size: 15x13x3.5cm, 0.133kgs
measuring range: 0-200cm

Packing: 1pc/box, 200pcs/ctn, 71.5x33x20cm, 18kgs



SW-G09
Tap measure

SW-G11

Syringe pen



SW-G39

Tongue depressor holder

MEDICAL GIFTS

NURSE WATCH

SW-G05A



SW-G05DW



SW-G05F



SW-G05JR



SW-G05KA



SW-G05KE



SW-G05KA



MEDICAL GIFTS

INSTRUMENTS BAG



SW-IB05

Instrument bag
these instruments is only samples packing for customer .
We can fitting these instruments following customer request.

SW-IB07



SW-IB06



SW-IB08



SW-IB09



MEDICAL GIFTS

PILL BOX



SW-PB28

7 room pill box



SW-PB32

7 room pill box



SW-PB34

28 ROOM PILL BOX, WITH CYLINDER DESIGN



SW-PB30

28 room pill box



SW-PB31

pill cut



SW-PB35

28 ROOM PILL BOX WITL LEATHER BAG



SW-PB33

7 room pill box