



HANGZHOU WEHERE MEDICAL TECHNOLOGY CO., LTD.

#701, Block C, Wanda Plaza, Gongshu, Hangzhou, 310015, China

Tel.: +86-571-86822701 Cell.: +8613905717534

E-mail: moneyxd@weheremed.com

www.weheremedical.com

WEHERE
MEDICAL
we here, we care



ANEROID SPHYGMOMANOMETER

PALM TYPE SPHYGMOMANOMETER

ELECTRONIC BLOOD PRESSURE MONITOR

STETHOSCOPE

HAMMER

PENLIGHT

PARTS

KS-1222

ECONOMY TYPE

Technical Features			
Aluminium alloy case			
Standard nylon cuff			
Standard Latex Free (PVC) bladder and bulb			
Dual use air-release valve with spring			
Taper end valve			
Vinyl zipper bag			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
100*70*195	520*400*370	50	20.5/19.5
Additional adult cuff color options			
■ Black ■ Gray ■ Blue			



KS-1621

STANDARD TYPE

Technical Features			
Aluminium alloy case			
Standard nylon cuff with D-ring			
Standard Latex Free (PVC) bladder and bulb			
Dual use air-release valve with spring			
Taper end valve			
Vinyl zipper bag			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
100*70*195	520*400*370	50	21/20
Additional adult cuff color options			
■ Black ■ Gray ■ Blue			



KS-1622

DELUXE TYPE

Technical Features			
Aluminium alloy case, red needle			
Standard nylon cuff			
Standard Latex Free (PVC) bladder and bulb			
Dual use air-release valve with spring			
Taper end valve			
Vinyl zipper bag			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
100*70*195	520*400*370	50	21/20
Additional adult cuff color options			
■ Black ■ Gray ■ Blue			



KS-1623

EURO TYPE

Technical Features			
Aluminium alloy case , PIN 80-160mmHg with red marking			
Adult Velcro blue calibrated cotton cuff			
Standard Latex Free (PVC) bladder and bulb			
Dual use air-release valve with spring			
Taper end valve			
Vinyl zipper bag			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
100*70*195	520*400*370	50	21/20
Additional adult cuff color options			
■ Black ■ Gray ■ Blue			



KS-3101

PREMIUM TYPE

Technical Features			
Nickel plated case			
Adult Velcro deluxe homecare nylon cuff with D-ring			
Standard Latex Free(PVC)bladder and bulb			
Dual use air-release valve with spring			
Vinyl zipper bag			
Included:adjustable key			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
100*70*195	520*400*370	50	21/20
Additional adult cuff color options			
■ Black ■ Gray ■ Blue			



KS-2205

DIY TYPE

Technical Features			
Aluminium alloy case			
Cotton cuff			
Latex bladder and bulb			
Dual use air-release valve with spring			
Taper end valve			
Vinyl zipper bag			
Double stethoscope			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
100*70*195	520*210*370	25	13.5/12.5
Additional adult cuff color options			
■ Black ■ Gray ■ Blue			



KS-2201

STANDARD SPHYGMOMANOMETER KIT

Technical Features			
Aluminium alloy case			
Standard nylon cuff			
Standard Latex Free (PVC) bladder and bulb			
Dual use air-release valve with spring			
Taper end valve			
Vinyl zipper bag			
Single stethoscope			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
100*70*195	520*210*370	25	11.5/10.5
Additional adult cuff color options			
Black Gray Blue			



KS-2021

HOME USE SPHYGMOMANOMETER KIT

Technical Features			
Aluminium alloy case			
Nylon cuff with D-ring , the single head stethoscope is fixed to cuff			
Standard Latex Free (PVC) bladder and bulb			
Dual use air-release valve with spring			
Taper end valve			
Vinyl zipper bag			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
100*70*195	520*210*370	25	11/10
Additional adult cuff color options			
Black Gray Blue			



KS-2203

SPHYGMOMANOMETER KIT WITH RAPPAORT STETHOSCOPE

Technical Features			
Aluminium alloy case			
Nylon cuff			
Standard Latex Free (PVC) bladder and bulb			
Dual use air-release valve with spring			
Taper end valve			
Shinning nylon carrying bag			
With rappaport stethoscope			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
140*70*230	715*245*370	25	18.5/17.5
Additional adult cuff color options			
Black Gray Blue			



KS-2201A

CHILD KIT TYPE

Technical Features			
Aluminium alloy case			
Classic cartoon new born cotton cuff + Classic cartoon child cotton cuff			
2-tube new born size pvc bladder+2-tube child size pvc bladder			
Standard pvc bulb			
Dual use air-release valve with spring			
Taper end valve			
Cartoon zipper bag			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
100*70*195	520*210*370	25	13/12



KS-2202

SPHYGMOMANOMETER KIT WITH SINGLE HEAD STETHOSCOPE

Technical Features			
Color aluminium alloy case			
Nylon color cuff			
Color Standard Latex Free (PVC) bladder and bulb			
Dual use air-release valve with spring			
Taper end valve			
Color Vinyl zipper bag			
With single stethoscope			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
100*70*195	520*210*370	25	11.5/10.5
Additional adult cuff color options			
Black Gray Blue			



KS-2203A

PRO'S COMBO CLINICIAN

Technical Features			
Aluminium alloy case			
Color-matched nylon cuff			
Standard Latex Free (PVC) bladder and bulb			
Dual use air-release valve with spring			
Taper end valve			
Leather zipper bag			
With rappaport stethoscope			
Nurse untily scissors,penlight,tape measure,and name tag			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
140*70*230	715*245*375	25	20.5/19.5
Additional adult cuff color options			
Black Gray Blue			



KS-3103

Technical Features			
Chrome plated metal case,The dial' s diameter is 48 mm			
Standard cotton cuff			
2-tube adult latex bladder			
Stainless steel handle with chrome plated			
PVC bulb with end valve			
Green nylon bag			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
100*70*195	520*210*370	25	17/16



KS-3202

Technical Features			
ABS plastic upper case			
The dial' s diameter is 60.2mm			
Stainless steel handle			
Cotton cuff			
1-tube pvc bladder			
Standard pvc blub			
Taper end valve			
Large vinyl zipper bag			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
100*70*195	520*210*370	25	12.5/11.5



KS-3102

Technical Features			
Chrome-plated housing			
The dial' s diameter is 48mm			
Stainless steel handle			
Nylon cuff			
2-tubes pvc bladder			
Standard latex blub			
Taper end valve			
Black nylon bag			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
100*70*195	520*210*370	25	13/12



2 tubes cuffs

KS-3301

Technical Features			
Abs housing			
The dial' s diameter is 48mm			
Stainless steel handle			
Nylon cuff			
1-tube pvc bladder			
Standard pvc blub			
Taper end valve			
Large vinyl zipper bag			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
100*70*195	520*210*370	25	12/11



KS-3302

Technical Features			
Nickel plated housing			
The dial' s diameter is 59.5mm			
Zinc-alloy handle with nickel plated			
Classic nylon cuff			
1-tube pvc bladder			
Ribbed pvc blub			
Taper end valve			
Large vinyl zipper bag			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
100*70*195	520*210*370	25	18/17



KS-3401

Technical Features			
Abs plastic housing , 0-300mmHg			
The dial' s diameter is 59.5mm			
Standard stainless steel handle			
Nylon cuff with D-ring			
2-tubes pvc bladder			
Standard pvc blub			
Taper end valve			
Large vinyl zipper bag			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
100*70*195	520*210*370	25	14/13



KS-3402

Technical Features			
Palm type sphygmomanometer			
Classic nylon cuff for infant/child/adult/large adult/thigh			
1-tube pvc bladder for infant/child/adult/large adult/thigh			
Standard pvc bulb			
Taper end valve			
Big carrying bag			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
460*90*255	480*470*280	5	7.6/6.6



KS-3203

Technical Features			
Square shaped dial 145mm*145mm			
Gradient adjustable for easy reading			
Standard nylon cuff			
2-tube adult pvc bladder			
Standard pvc blub			
Taper end valve			
Spiral pvc tube(1.8m)			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
185*160*180	850*390*380	20	21/20



KS-3203A

Technical Features			
Square shaped dial 145mm*145mm			
360°swivel for easy reading			
Gradient adjustable for easy reading			
5-wheel with brake			
Adjustable height			
Metal basket			
Nylon cuff			
2-tube adult pvc bladder+pvc bulb			
Spiral pvc tube(1.8m)			
Plastic connector			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
/	505*480*185	1	3.2/2.2



Electronic blood pressure monitor

KS-0401

Technical Features			
LCD display			
Measurement method:oscillometric method			
Rated voltage:DC 6V(4AAA batteries)/D.C.3.7V			
Power supply:AAA battery/lithium battery			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
100*80*160	530*430*345	50	17.5/16.5



KS-0402

Technical Features			
LED display			
Measurement method:oscillometric method			
Rated voltage:DC 6V(4AAA batteries)/D.C.3.7V			
Power supply:AAA battery/lithium battery			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
100*80*160	530*430*345	50	21.5/20.5



KS-0403

Technical Features			
LD display			
Measurement method:oscillometric method			
Rated voltage:DC 3.7V			
Power supply:lithium battery			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
100*80*80	520*440*350	100	21/20



KS-2035[Adult] KS-2036[Child]

Technical Features			
Single head chestpiece for adult.aluminum chestpiece,the ϕ 43mm membrane			
Applicable for ausscultating of heart and lungs			
Single head chestpiece for child.aluminum chestpiece,the ϕ 33mm membrane			
Applicable for ausscultating of heart and lungs			
The small bell-like chestpiece is used for auscultating slight and low frequency sounds			
Standard single pvc tube,color is optional			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
100*25*200	520*380*420	100	11/10



KS-2035A

Technical Features			
Single head chestpiece for adult.aluminum color chestpiece,the ϕ 43mm membrane			
Color aluminum binaural			
Applicable for ausscultating of heart and lungs			
The small bell-like chestpiece is used for auscultating slight and low frequency sounds			
Standard single pvc tube,color is optional			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
100*25*200	520*380*420	100	11/10



KS-2032A

Technical Features			
Abs chestpiece			
PET binaural			
Pvc tube			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
Self sealing bag	440*400*310	100	6/5



KS-2045[Adult] KS-2046[Child]

Technical Features			
Dual head chestpiece for adult.aluminum chestpiece,the ϕ 43mm membrane			
Applicable for ausscultating of heart and lungs			
Dual head chestpiece for child.aluminum chestpiece,the ϕ 33mm membrane			
Applicable for ausscultating of heart and lungs			
The small bell-like chestpiece is used for auscultating slight and low frequency sounds			
Standard single pvc tube,color is optional			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
100*25*200	520*380*420	100	12/11



KS-2045A

Technical Features			
Dual head chestpiece for adult.aluminum color chestpiece,the ϕ 43mm membrane			
Color aluminum binaural			
Applicable for ausscultating of heart and lungs			
The small bell-like chestpiece is used for auscultating slight and low frequency sounds			
Standard single pvc tube,color is optional			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
100*25*200	520*380*420	100	12/11



KS-2045B

Technical Features			
Dual head chestpiece for adult.aluminum color chestpiece,with color-matched non-chill rim			
Aluminum binaural			
Applicable for ausscultating of heart and lungs			
The small bell-like chestpiece is used for auscultating slight and low frequency sounds			
Standard single pvc tube,color is optional			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
100*25*200	520*380*420	100	12/11



KS-2013

Technical Features			
Silver dual head chestpiece for neonate and baby			
Applicable for auscultating of the hearts and lungs of neonate and baby			
Zinc alloy chestpiece,good sound conduction			
Considerately dual anti-chill ring design			
Standard single pvc tube,color is optional			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
100*25*200	520*380*420	100	14.5/13.5



KS-2015

Technical Features					
Silver dual-membrane chestpiece:the big membrane with a diameter of 45.5mm is sensitive for sound collecting,the small membrane is applicable for child					
Zinc alloy chestpiece,good sound conduction					
Standard dual pvc tube					
Attached transparent multi-pocket bag included					
Package					
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn		
296*130*40	670*310*430	50	19/18		
Additional adult cuff color options					
■ Black	■ Gray	■ Blue	■ Pink	■ Yellow	■ Green
■ Red	■ Purple	■ Orange	■ Wine red		



KS-2019

- Black Stealth chestpiece and binaural



KS-2018

- Gold-plated chestpiece and binaural

KS-2038

Technical Features			
Clock function,special design			
Single chestpiece:the membrane with a diameter of 45.5mm			
Zinc alloy chestpiece,good sound conduction			
Standard dual pvc tube			
Attached transparent multi-pocket bag included			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
296*130*40	670*310*430	50	20/19



KS-2063

Technical Features			
Single head stethoscope fpr teaching use with two binaural			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
296*130*40	670*310*430	50	12/11



KS-2027(Child) KS-2026(Adult)

Technical Features			
Silver dual head chestpiece for adult			
The small bell-like chestpiece is used for auscultating slight and low frequency sounds			
Considerately dual anti-chill ring design			
Stainless steel binaural			
Stainless steel chestpiece ,excellent sound conduction			
Includes : A pair of metal ear tips and a replaceable membrane			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
165*50*320	840*335*270	25	13/12



KS-2026B

Technical Features			
Silver dual head chestpiece for adult,Black			
The small bell-like chestpiece is used for auscultating slight and low frequency sounds			
Considerately dual anti-chill ring design			
Stainless steel black binaural			
Stainless steel black chestpiece ,excellent sound conduction			
Includes : A pair of metal ear tips and a replaceable membrane			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
165*50*320	840*335*270	25	13/12



KS-2026BC

Technical Features			
Bronze stainless steel chestpiece			
The small bell-like chestpiece is used for auscultating slight and low frequency sounds			
Considerately dual anti-chill ring design			
Bronze stainless steel binaural			
Customized patterned tubes are available			
Includes:A pair of metal ear tips and a replaceable membrane			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
165*50*320	840*335*270	25	13.2/12.2



KS-2026BT

Technical Features			
Colorful stainless steel chestpiece			
The small bell-like chestpiece is used for auscultating slight and low frequency sounds			
Considerately dual anti-chill ring design			
Colorful stainless steel binaural			
Customized patterned tubes are available			
Includes:A pair of metal ear tips and a replaceable membrane			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
165*50*320	840*335*270	25	13.2/12.2



KS-2026BG

Technical Features			
Gold-plated stainless steel binaural			
The small bell-like chestpiece is used for auscultating slight and low frequency sounds			
Considerately dual anti-chill ring design			
Gold-plated stainless steel binaural			
Includes:A pair of metal ear tips and a replaceable membrane			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
165*50*320	840*335*270	25	13/12



KS-2036A

Technical Features			
Stainless steel chestpiece			
The small bell-like chestpiece is used for auscultating slight and low frequency sounds			
Considerately dual anti-chill ring design			
Stainless steel binaural			
Includes:A pair of metal ear tips and a replaceable membrane			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
165*50*320	840*335*270	25	13/12



Float diaphragm

KS-2036B

Technical Features			
Black stainless steel chestpiece			
The small bell-like chestpiece is used for auscultating slight and low frequency sounds			
Considerately dual anti-chill ring design			
Black stainless steel binaural			
Includes:A pair of metal ear tips and a replaceable membrane			
Package			
Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
165*50*320	840*335*270	25	13/12



KS-2025A

Technical Features

Silver dual head chestpiece for cardiology diagnosis
The small ,deep,cone-shaped chestpiece is used for auscultating of cacophony, as well as slight and low frequency sounds
Stainless steel chestpiece,excellent sound conduction
Considerately dual anti-chill ring design
"A" shaped tube ensures the best sound conduction
Connecting spring is embedded ,ear tips are specially designed for suiting the ear structure
Thick tube with excellent resistance to ambient interference
A replaceable stainless steel chestpiece
Includes:A pair of ear tips,and a replaceable float diaphragm

Package

Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
165*50*320	840*335*270	25	14/13



KS-2025

Technical Features

Silver dual head chestpiece for cardiology diagnosis
The small ,deep,cone-shaped chestpiece is used for auscultating of cacophony, as well as slight and low frequency sounds
Zinc alloy chestpiece , excellent sound conduction
Considerately dual anti-chill ring design
Connecting spring is embedded ,ear tips are specially designed for suiting the ear structure
Thick tube with excellent resistance to ambient interference
Includes:A pair of ear tips,and a replaceable float diaphragm

Package

Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
250*38*132	675*265*420	50	14.8/13.8



KS-2055

Technical Features

Silver single head chestpiece,ergonomic design for comfortable use
The diameter of chestpiece is 43mm,and it can catch both low frequency and high frequency sounds by pressing the back
Zinc alloy chestpiece with nickel-plated silver finish,excellent sound conduction, little acoustoc attenuation
Connecting spring is embedded,ear tips are specially designed to fit the ear structure
Considerately anti-chill ring design

Package

Box size(mm)	Carton size(mm)	Pcs/Ctn	Kgs/Ctn
250*38*132	675*265*420	50	15/14



Hammer

KS-0201

Neurological Hammer



Weight	Length	Width
56 (g)	188 (mm)	21 (mm)
Carton size	PCS/CTN	KGS/CTN
570*210*230	100	8.5
Packing bill		
Blister		
Annotation :		
Zinc alloy handle , PVC head , color is optional		

KS-0202

Neurological Hammer

Weight	Length	Width
56 (g)	188 (mm)	21 (mm)
Carton size	PCS/CTN	KGS/CTN
635*240*205	100	10
Packing bill		
Blister		
Annotation :		
Zinc alloy handle , PVC head , Diagnostic Pin can be screwed off conveniently in the head , while diagnostic brush can be taken out in the end.		



KS-0203


Neurological Hammer



Weight	Length	Width
60g)	190 (mm)	55 (mm)
Carton size	PCS/CTN	KGS/CTN
520*200*220	100	9
Packing bill		Box size
Blister		190*55*20mm
Annotation :		
Zinc alloy handle , PVC head , color is optional		

KS-0204

Medical gift



Carton size	PCS/CTN	KGS/CTN
570*365*300mm	40	23
Packing bill		Box size
Plastic box		345*130*30mm
Annotation :		
Diagnostic kit		

KS-0204

Stainless Steel Penlight

Weight	Length	Width
16 (g)	136 (mm)	12 (mm)
Carton size	PCS/CTN	KGS/CTN
380*330*300	400	9.8
Packing bill		Box size
Blister		145*28*16mm
Annotation :		
2×AAA batteries used		



KS-0302

Stainless Steel Penlight



Weight	Length	Width
14.3 (g)	136 (mm)	12 (mm)
Carton size	PCS/CTN	KGS/CTN
510*370*240mm	250	16
Packing bill		Box size
Gift box		175*160*210mm
Annotation :		
2×AAA batteries used , stainless steel shell		

Chestpiece



Gauge		
01		
specification	Model No	Surface Treatment
Special Aluminum alloy case	01	Powder-coating finished
Carton size	Pcs/ctn	Kgs/ctn
42*29.5*27.5cm	100	13/12

Bulb		
02pvc ribbed bulb	03standard latex bulb	04standard pvc bulb

Case (logo printing is available)			
05vinyl zipper case	06nylon bag	07nylon carrying case	08black drawstring bag

Bladder



09 latex bladder with 1-tube



10 latex bladder with 2-tube



11 pvc bladder with 1-tube



12 pvc bladder with 2-tube

Cuff



13 standard nylon cuff



14 standard cotton cuff



15 D-ring cuff



16 cuff with head attached

Binaural



17 grooved aluminum binaural



18 whorl stainless steel binaural

Tube



19 Colored PVC Tube



20 PVC Tube for Rappaport Type



21 PVC Transparent Tube