

(0086571-86822701



Add.: # 701, Building C, Wanda Plaza, No. 666,
Hangxing Road, Gongshu District,
Hangzhou, 310015, Zhejiang Province, China

Web: www.weheremedical.com

Wehere

Hangzhou Wehere Medical Technology Co., Ltd.

Umbilical Catheter



Umbilical Catheter

Catheter in radiopaque and transparent Polyurethane, can be used for umbilical venous catheterisation(parenteral nutrition, drug administration, exchange transfusion, blood sampling...) or for umbilical arterial catheterisation (blood gas analysis, blood pressure measurement...)

Benefits:

- Flexible, thermos-sensitive material with open and rounded distal tip to reduce vascular trauma
- DEHP-free
- Polyurethane has good haemocompatibility: reduces risk of thrombosis



Single Lumen				
Catalogue	Lumen	TypeLumen	Length(cm)	3port
UC1-20-30	Single Lumen	20Ga	30	×
UC1-20-30-W	Single Lumen	20Ga	30	√
UC1-18-30	Single Lumen	18Ga	30	×
UC1-18-30-W	Single Lumen	18Ga	30	√
UC1-16-30	Single Lumen	16Ga	30	×
UC1-16-30-W	Single Lumen	16Ga	30	√
Double Lumen				
Catalogue	Lumen	TypeLumen	Length(cm)	3port
UC2-4-30	Double Lumen	4F	30	×
UC2-4-30-W	Double Lumen	4F	30	√
UC2-5-30	Double Lumen	5F	30	×
UC2-5-30-W	Double Lumen	5F	30	√
Triple Lumen				
Catalogue	Lumen	TypeLumen	Length(cm)	3port
UC3-4-30	Triple Lumen	4F	30	×
UC3-4-30-W	Triple Lumen	4F	30	√
UC3-5-30	Triple Lumen	5F	30	×
UC3-5-30-W	Triple Lumen	5F	30	√