

(0086571-86822701


Add.: # 701, Building C, Wanda Plaza, No. 666,
Hangxing Road, Gongshu District,
Hangzhou, 310015, Zhejiang Province, China

Web: www.weheremedical.com

Wehere

Hangzhou Wehere Medical Technology Co., Ltd.

Arterial Catheter



Arterial Catheter



Product Introduction

An arterial catheter (also called an arterial line or A-line) is a thin,flexible tube inserted into an artery to monitor continuous bloodpressure, obtain frequent blood samples, or assess hemodynamicstatus in critically ill patients.

Product Features

- Real-time blood pressure monitoring (more accurate than non-invasive cuffs).
- Direct arterial blood sampling (forABGs,labs).
- Common insertion sites: Radial artery (most common), femoral, brachial, or dorsalis pedis.
- **Uses:** ICU, major surgery, shock, hemodynamic instability.

Product Specifications

Product No	Size
681001	18G
681002	20G
681003	22G
681004	24G