



HANGZHOU WEHERE MEDICAL TECHNOLOGY CO., LTD.

#701, Block C, Wanda Plaza, Gongshu, Hangzhou, 310015, China

Tel.: +86-571-86822701 Cell.: +8613905717534

E-mail: moneyxd@weheremed.com

www.weheremedical.com

WEHERE
MEDICAL
we here, we care




两件式产品剖析图

TWO-PIECE PRODUCT ANALYSIS PROFILE


离型纸 (Release paper)

吸附性强，易掀开
Great adhesive, easy to peel off




子母安全扣件 (Flange)

轻松分离，安全牢靠
搭配腰带，密封更好
Easy to be removed separately, safe and secure, leak proof if combines with belt



活性炭片 (Carbon filter)

吸附性强，有效祛除异味
Strong absorption, neutralize odor



棕色袋体 (Brown bag)

阻隔性强，隐蔽性好，避免尴尬
High barrier, good privacy, avoid embarrassment



半透明水胶体
Translucent hydro-colloid

采用亲肤型水胶体材质，持久亲肤，不易过敏，粘贴无残留
Skin friendly hydro-colloid material, long lasting and skin friendly, no stimulation, no residual on the skin



医用无纺布
Medical non-woven fabric

柔软舒适透气，高阻隔静音膜，降低摩擦声，避免尴尬
Soft, comfortable and ventilate, high barrier and silent film, reduce friction noise, avoid embarrassment



01

国内系列 (Domestic Range)

一件式肛肠造口袋

One-piece Colostomy Bag

产品描述 (Product Description)

- 本产品为回肠、结肠、直肠造口的护理产品。
- 产品接触完好皮肤和肠内腔。
- 非无菌提供。
- This is a medical product for ileum, colon and rectum.
- Product comes into contact with intact skin and intestinal lumen.
- Supplied non-sterile.

产品简介 (Product Introduction)

项目	介绍
组成结构	主要由剪孔离型纸（膜）、水胶体底盘、高阻隔膜袋体、无纺布、活性炭和封口器组成。
产品特点	操作使用方便，底盘黏性强，对皮肤友好，不易过敏；袋体内置活性炭片，有效减少异味。
Composition Parts	It is made of pre-cut PET paper (film), hydrocolloid base plate, EVOH film, non-woven fabric, carbon filter and drainable end.
Product Features	Easy to operate, great adhesive base plate, skin friendly, rarely get allergy; inserted carbon filter, avoid embarrassment effectively.

型号(Model): 1106
剪孔(Cut to fit): 15~65mm

封口方式:扎带封口
底盘直径:108mm
无纺布颜色:白色
袋体颜色:黄色+透明

Medical twist tie
Skin barrier diameter: 108mm
White non-woven fabric
Yellow and transparent bag



型号(Model): 1102
剪孔(Cut to fit): 15~55mm

封口方式:维克罗搭扣
底盘直径:108mm
无纺布颜色:白色
袋体颜色:双面茶色

Velcro
Skin barrier diameter: 108mm
White non-woven fabric
Double sides brown bag



型号(Model): 1105
剪孔(Cut to fit): 15~85mm

封口方式:扎带封口
底盘直径:125mm
无纺布颜色:白色
袋体颜色:双面茶色

Medical twist tie
Skin barrier diameter: 125mm
White non-woven fabric
Double sides brown bag



型号(Model): 1169
剪孔(Cut to fit): 15~70mm

封口方式:扎带封口
底盘直径:108mm
无纺布颜色:白色
袋体颜色:黄色+透明

Medical twist tie
Skin barrier diameter: 108mm
White non-woven fabric
Yellow and transparent bag



型号(Model): 1103
剪孔(Cut to fit): 15~60mm

封口方式:夹子封口
底盘直径:108mm
无纺布颜色:白色
袋体颜色:双面茶色

Clamp
Skin barrier diameter: 108mm
White non-woven fabric
Double sides brown bag



型号(Model): 1109
剪孔(Cut to fit): 15~60mm

封口方式:夹子封口
底盘直径:108mm
无纺布颜色:白色
袋体颜色:双面茶色

Clamp
Skin barrier diameter: 108mm
White non-woven fabric
Double sides brown bag



01-1

国内系列 (Domestic Range)

一件式肛肠造口袋

One-piece Colostomy Bag

产品描述 (Product Description)

- 本产品为回肠、结肠、直肠造口的护理产品。
- 产品接触完好皮肤和肠内腔。
- 非无菌提供。
- This is a medical product for ileum, colon and rectum.
- Product comes into contact with intact skin and intestinal lumen.
- Supplied non-sterile.

产品简介 (Product Introduction)

项目	介绍
组成结构	主要由剪孔离型纸（膜）、水胶体底盘、高阻隔膜袋体、无纺布、活性炭和封口器组成。
产品特点	操作使用方便，底盘黏性强，对皮肤友好，不易过敏；袋体内置活性炭片，有效减少异味。
Composition Parts	It is made of pre-cut PET paper (film), hydrocolloid base plate, EVOH film, non-woven fabric, carbon filter and drainable end.
Product Features	Easy to operate, great adhesive base plate, skin friendly, rarely get allergy; inserted carbon filter, avoid embarrassment effectively.

型号(Model): 1107
剪孔(Cut to fit): 15~65mm

封口方式:魔术贴
底盘直径:108mm
无纺布颜色:白色
袋体颜色:黄色+透明

Hook and loop
Skin barrier diameter: 108mm
White non-woven fabric
Yellow bag and transparent bag





型号(Model): 1338
剪孔(Cut to fit): 15~38mm

封口方式:维克罗搭扣
凸面底盘直径:108mm
无纺布颜色:白色
袋体颜色:双面透明

Velcro
Convex skin barrier diameter: 108mm
White non-woven fabric
Double sides transparent bag





型号(Model): 1104
剪孔(Cut to fit): 15~48mm

封口方式:扎带封口
底盘直径:87mm
无纺布颜色:白色
袋体颜色:双面透明

Medical twist tie
Skin barrier diameter: 87mm
White non-woven fabric
Double sides transparent bag





型号(Model): 1209
剪孔(Cut to fit): 15~60mm

封口方式:夹子封口
底盘直径:99mm
底盘类型:双层底盘
袋体颜色:黄色+透明

Clamp
Skin barrier diameter: 99mm
With spunlace extender
Yellow bag and transparent bag





型号(Model): 1261
剪孔(Cut to fit): 15~57mm

封口方式:夹子封口
底盘直径:87mm
底盘类型:双层底盘
袋体颜色:双面茶色

Clamp
Skin barrier diameter: 87mm
With spunlace extender
Double sides brown bag





型号(Model): 1130
剪孔(Cut to fit): 15~70mm

封口方式:魔术贴
底盘直径:118mm
灰色防水涤纶布
袋体颜色:双面灰色

Hook and loop
Skin barrier diameter: 118mm
Gray Waterproof Fabric
Gray bag





01-2

国内系列 (Domestic Range)

一件式造口袋 One-piece Ostomy Bag

产品描述 (Product Description)

- 本产品为回肠、结肠、直肠造口的护理产品。
- 产品接触完好皮肤和肠内腔。
- 非无菌提供。
- This is a medical product for ileum, colon and rectum.
- Product comes into contact with intact skin and intestinal lumen.
- Supplied non-sterile.

产品简介 (Product Introduction)

项目	介绍
组成结构	主要由剪孔离型纸（膜）、水胶体底盘、高阻隔膜袋体、无纺布、活性炭和封口器组成。
产品特点	操作使用方便，底盘黏性强，对皮肤友好，不易过敏；袋体内置活性炭片，有效减少异味。
Composition Parts	It is made of pre-cut PET paper (film), hydrocolloid base plate, EVOH film, non-woven fabric, carbon filter and drainable end.
Product Features	Easy to operate, great adhesive base plate, skin friendly, rarely get allergy; inserted carbon filter, avoid embarrassment effectively.

型号(Model): 3109
剪孔(Cut to fit): 15~50mm

袋体类型:泌尿造口袋
底盘直径:108mm
无纺布颜色:白色
防逆流

One-piece Urostomy bag
Skin barrier diameter: 108mm
White non-woven fabric
Anti-backflow



型号(Model): 3105
剪孔(Cut to fit): 15~45mm

袋体类型:泌尿造口袋
底盘直径:108mm
无纺布颜色:白色
防逆流

One-piece Urostomy bag
Skin barrier diameter: 108mm
White non-woven fabric
Anti-backflow



型号(Model): 3205
剪孔(Cut to fit): 15~48mm

袋体类型:泌尿造口袋
底盘直径:87mm
底盘类型:双层底盘
防逆流

One-piece Urostomy bag
Skin barrier diameter: 87mm
With spunlace extender
Anti-backflow



型号(Model): 1301
规格(Volume): 1500ml

袋体类型:集尿袋
最大容量:1500ml
防逆流型

Urine collection bag
Max Volume: 1500ml
Anti-backflow



型号(Model): 2101
剪孔(Cut to fit): 15~50mm

袋体类型:肛肠闭口袋
底盘直径:92mm
无纺布颜色:白色
袋体颜色:双面茶色

Closed bag
Skin barrier diameter: 92mm
White non-woven fabric
Double sides brown bag



型号(Model): 2102
剪孔(Cut to fit): 15~35mm

造口袋类型:闭口袋
底盘直径:87mm
无纺布颜色:白色
袋体颜色:黄色+透明

Closed bag
Skin barrier diameter: 87mm
White non-woven fabric
Yellow and transparent bag



01-3



国内系列

(Domestic Range)

二件式造口袋

Two-piece Ostomy Bag

产品描述

(Product Description)

- 本产品为回肠、结肠、直肠造口的护理产品。

- 产品接触完好皮肤和肠内腔。

- 非无菌提供。

- This is a medical product for ileum, colon and rectum.

- Product comes into contact with intact skin and intestinal lumen.

- Supplied non-sterile.



产品简介

(Product Introduction)

项目	介绍
组成结构	底盘(7250)由离型膜、水胶体、扣件系统组成；袋体(4256)主要由高阻隔膜袋体、无纺布、活性碳片、扣件系统和封口器组成。
产品特点	灰色袋体隐蔽性好，底盘黏性强，对皮肤友好，锁扣式设计，密封性更好，使用更安全。

Item	Introduction
Components	Skin barrier (7250) is composed of PET film, hydrocolloid and flange system; The bag (4256) is mainly composed of EVOH film, non-woven fabric, carbon filter, flange and drainable end.
Product Features	The gray bag has good concealment, the skin barrier is with strong adhesive and skin-friendly, the lock-type design provides better sealing and even safer.

产品简介

(Product Introduction)

项目	介绍
组成结构	造口底盘(7338&7301)由离型纸、水胶体、扣件系统组成；袋体(4209&4301)主要由高阻隔膜袋体、无纺布、活性碳片、扣件系统和封口器组成。
产品特点	操作使用方便，底盘黏性强，对皮肤友好，不易过敏；袋内置活性炭片，有效减少异味。

Item	Introduction
Components	The stoma skin barrier (7338&7301) is composed of paper, hydrocolloid and flange system; The bag (4201&4301) is mainly composed of EVOH film, non-woven fabric, carbon filter, flange and drainable end .
Product Features	Easy to operate and use, strong adhesiveness, skin friendly, and not easy to get allergic; built-in carbon filter, which can effectively reduce odor.

02

医院系列 (Hospital Range)

一件式造口袋 One-piece Ostomy Bag

产品描述 (Product Description)

- 本产品为回肠、结肠、直肠造口的护理产品。
- 产品接触完好皮肤和肠内腔。
- 非无菌提供。
- This is a medical product for ileum, colon and rectum.
- Product comes into contact with intact skin and intestinal lumen.
- Supplied non-sterile.

产品简介 (Product Introduction)

项目	介绍
组成结构	主要由剪孔离型纸（膜）、水胶体底盘、高阻隔膜袋体、无纺布、活性炭和封口器组成。
产品特点	操作使用方便，底盘黏性强，对皮肤友好，不易过敏；袋体内置活性炭片，有效减少异味。
Composition Parts	It is made of pre-cut PET paper (film), hydrocolloid base plate, EVOH film, non-woven fabric, carbon filter and drainable end.
Product Features	Easy to operate, great adhesive base plate, skin friendly, rarely get allergy; inserted carbon filter, avoid embarrassment effectively.

型号(Model): 3171
剪孔(Cut to fit): 15~60mm

袋体类型:泌尿造口袋
底盘直径:99mm
无纺布颜色:白色
防逆流


One-piece Urostomy bag
Skin barrier diameter: 99mm
White non-woven fabric
Anti-backflow



型号(Model): 1173
剪孔(Cut to fit): 15~60mm

封口方式:夹子
底盘直径:108mm
无纺布颜色:白色
背开窗


Clamp
Skin barrier diameter: 108mm
White non-woven fabric
Back opening window



型号(Model): 1174
剪孔(Cut to fit): 15~35mm

封口方式:扎带封口
底盘直径:87mm
无纺布颜色:背开窗设计
背开窗

Medical twist tie
Skin barrier diameter: 87mm
White non-woven fabric
Back opening window



型号(Model): 3172
剪孔(Cut to fit): 15~45mm

袋体类型:泌尿造口袋
底盘直径:108mm
无纺布颜色:白色
防逆流


One-piece Urostomy bag
Skin barrier diameter: 108mm
White non-woven fabric
Anti-backflow



型号(Model): 1175
剪孔(Cut to fit): 15~85mm

封口方式:维克罗搭扣
底盘直径:125mm
无纺布颜色:白色
背开窗

Velcro
Skin barrier diameter: 125mm
White non-woven fabric
Back opening window



型号(Model): 2172
剪孔(Cut to fit): 15~35mm

造口袋类型:闭口袋
底盘直径:87mm
无纺布颜色:白色
袋体颜色:透明

Closed bag
Skin barrier diameter: 87mm
White non-woven fabric
Transparent bag



02-1

医院系列

(Hospital Range)

二件式造口底盘

Two-piece Skin Barrier

- 产品描述 (Product Description)
- 本产品为回肠、结肠、直肠造口的护理用品。
 - 产品接触完好皮肤和肠内腔。
 - 非无菌提供。
 - This is a medical product for ileum, colon and rectum.
 - Product comes into contact with intact skin and intestinal lumen.
 - Supplied non-sterile.

型号(Model): 7331 (凸面)

剪孔(Cut to fit): 15~38mm

扣件(Flange): 45扣件



型号(Model): 7332 (凸面)

剪孔(Cut to fit): 15~45mm

扣件(Flange): 57扣件



型号(Model): 7271

剪孔(Cut to fit): 15~45mm

扣件(Flange): 45扣件



型号(Model): 7371

剪孔(Cut to fit): 15~57mm

扣件(Flange): 57扣件



型号(Model): 7471

剪孔(Cut to fit): 15~70mm

扣件(Flange): 70扣件



■ 产品简介 (Product Introduction)

项目	介绍	Item	Introduction
组成结构	二件式凸面底盘由进口水胶体材质、离型膜和扣件系统组成；	Components	The two-piece convex skin barrier is composed of imported hydrocolloid material, PET film and flange system;
产品特点	适用于造口周围皮肤凹陷者使用，独立包装，粘性好，对皮肤致敏率低，无刺激。扣环与袋体相对分离的设计，与袋体扣合和拆解时更方便。	Product Features	It is suitable for those with sunken skin around the stoma, individually packaged, strong adhesiveness, low skin sensitization rate, and no irritation. The design that the flange is accordingly separated from the bag is more convenient when it is fastened and disassembled with the bag.

■ 产品简介 (Product Introduction)

项目	介绍	Item	Introduction
组成结构	二件式底盘由进口水胶体材质、离型膜和扣件系统组成；	Components	The two-piece skin barrier is composed of imported hydrocolloid material, PET film and flange system;
产品特点	独立包装的二件式底盘，粘性好，对皮肤致敏率低，无刺激。扣环与袋体相对分离的设计，与袋体扣合和拆解时更方便。	Product Features	Individually packaged two-piece skin barrier, strong adhesiveness, low skin sensitization rate, and no irritation. The design that the flange is accordingly separated from the bag is more convenient when it is fastened and disassembled with the bag.

02-2



医院系列

(Hospital Range)

二件式造口袋

Two-piece Ostomy Bag

■ 产品描述 (Product Description)

- 本产品为回肠、结肠、直肠造口的护理用品。

- 产品接触完好皮肤和肠内腔。

- 非无菌提供。

- This is a medical product for ileum, colon and rectum.

- Product comes into contact with intact skin and intestinal lumen.

- Supplied non-sterile.



■ 产品简介 (Product Introduction)

项目	介绍	Item	Introduction
组成结构	袋体膜材为新型高分子双面透明薄膜，高阻隔材料，能有效隔绝异味。单面无纺布内衬，柔软亲肤，减少摩擦外衬全透明方便观察造口情况和尿液状况。	Components	The bag film is a new type of EVOH double sides transparent film, which can effectively isolate odors. The single side non-woven lining is soft and skin-friendly, reducing friction and the outer lining is fully transparent to facilitate observation of stoma and urine conditions.
产品特点	适用于尿路造口者使用，防逆流设计，有效减少尿液对造口和底盘的腐蚀。	Product Features	It is suitable for urostomy patients. The anti back flow design can effectively reduce the corrosion of urine to the stoma and the skin barrier.

■ 产品简介 (Product Introduction)

项目	介绍	Item	Introduction
组成结构	袋体膜材为新型高分子双面透明薄膜，高阻隔材料，能有效隔绝异味； 双面无纺布，背开窗设计，柔软亲肤减摩擦。	Components	It is made of pre-cut PET paper (film), hydrocolloid base plate, EVOH film, non-woven fabric, carbon filter and drainable end.
产品特点	开口袋可随时排放内容物并拆解下来清洗再使用，方便实惠。 闭口袋轻便小巧，方便携带，干净更换。	Product Features	The drainable ostomy bag can discharge the contents at any time and disassemble it for cleaning and reuse, which is convenient and affordable. The closed bags are light, portable, clean and convenient.

03

外贸系列 (International Range)

一件式肛肠造口袋
One-piece Colostomy Bag

产品描述 (Product Description)

- 本产品为回肠、结肠、直肠造口的护理产品。
- 产品接触完好皮肤和肠内腔。
- 非无菌提供。
- This is a medical product for ileum, colon and rectum.
- Product comes into contact with intact skin and intestinal lumen.
- Supplied non-sterile.

产品简介 (Product Introduction)

项目	介绍
组成结构	主要由剪孔离型纸（膜）、水胶体底盘、高阻隔膜袋体、无纺布、活性炭和封口器组成。
产品特点	操作使用方便，底盘黏性强，对皮肤友好，不易过敏；袋体内置活性炭片，有效减少异味。
Composition Parts	It is made of pre-cut PET paper (film), hydrocolloid base plate, EVOH film, non-woven fabric, carbon filter and drainable end.
Product Features	Easy to operate, great adhesive base plate, skin friendly, rarely get allergy; inserted carbon filter, avoid embarrassment effectively.

型号(Model): Y1166
剪孔(Cut to fit): 15~70mm

封口方式: 维克罗搭扣
底盘直径: 118mm
无纺布颜色: 棕色
袋体颜色: 茶色

Velcro
Skin barrier diameter: 118mm
Brown non-woven fabric
Brown bag



型号(Model): Y1216
剪孔(Cut to fit): 15~60mm

封口方式: 夹子
底盘直径: 99mm
无纺布颜色: 棕色
袋体颜色: 茶色

Clamp
Skin barrier diameter: 99mm
Brown non-woven fabric
Brown bag



型号(Model): Y1120
剪孔(Cut to fit): 15~75mm

封口方式: 维克罗搭扣
底盘直径: 118mm
无纺布颜色: 白色
袋体颜色: 茶色

Velcro
Skin barrier diameter: 118mm
White non-woven fabric
Brown bag



型号(Model): Y1188
剪孔(Cut to fit): 15~105mm

封口方式: 夹子
底盘直径: 140mm
无纺布颜色: 白色
袋体颜色: 透明

Clamp
Skin barrier diameter: 140mm
White non-woven fabric
Transparent bag



型号(Model): Y1189
剪孔(Cut to fit): 15~60mm

封口方式: 夹子
底盘直径: 99mm
袋体颜色: 透明

Clamp
Skin barrier diameter: 99mm
Transparent bag



型号(Model): 1154
剪孔(Cut to fit): 15~35mm

封口方式: 扎带封口
无纺布颜色: 白色
袋体颜色: 透明

Medical twist tie
White non-woven fabric
Transparent bag



03-1

外贸系列 (International Range)

一件式造口袋 One-piece Ostomy Bag

产品描述 (Product Description)

- 本产品为回肠、结肠、直肠造口的护理产品。
- 产品接触完好皮肤和肠内腔。
- 非无菌提供。
- This is a medical product for ileum, colon and rectum.
- Product comes into contact with intact skin and intestinal lumen.
- Supplied non-sterile.

产品简介 (Product Introduction)

项目	介绍
组成结构	主要由剪孔离型纸（膜）、水胶体底盘、高阻隔膜袋体、无纺布、活性炭和封口器组成。
产品特点	操作使用方便，底盘黏性强，对皮肤友好，不易过敏；袋体内置活性炭片，有效减少异味。
Composition Parts	It is made of pre-cut PET paper (film), hydrocolloid base plate, EVOH film, non-woven fabric, carbon filter and drainable end.
Product Features	Easy to operate, great adhesive base plate, skin friendly, rarely get allergy; inserted carbon filter, avoid embarrassment effectively.

型号(Model): Y3103
剪孔(Cut to fit): 15~60mm

袋体类型:泌尿造口袋
底盘直径:108mm
无纺布颜色:白色
防逆流

One-piece Urostomy bag
Skin barrier diameter: 108mm
White non-woven fabric
Anti-backflow



型号(Model): Y1105
剪孔(Cut to fit): 15~85mm

封口方式:夹子
底盘直径:125mm
无纺布颜色:白色
袋体颜色:双面茶色

Clamp
Skin barrier diameter: 125mm
White non-woven fabric
Brown bag



型号(Model): Y1118B
剪孔(Cut to fit): 10~75mm

封口方式:魔术贴
底盘直径:Max 110mm
无纺布颜色:棕色
袋体颜色:黄色+透明

Hook and loop
Skin barrier diameter: Max 110mm
Brown non-woven fabric
Yellow and transparent bag



型号(Model): Y3201
剪孔(Cut to fit): 15~48mm

袋体类型:泌尿造口袋
底盘直径:87mm
底盘类型:双层底盘
防逆流

One-piece Urostomy bag
Skin barrier diameter: 87mm
With spunlace extender
Anti-backflow



型号(Model): Y1203
剪孔(Cut to fit): 15~60mm

封口方式:维克罗搭扣
底盘直径:99mm
底盘类型:双层底盘
袋体颜色:茶色

Velcro
Skin barrier diameter: 99mm
With spunlace extender
Brown bag



型号(Model): Y1124
剪孔(Cut to fit): 15~70mm

封口方式:夹子
底盘直径:118mm
无纺布颜色:白色
袋体颜色:黄色+透明

Clamp
Skin barrier diameter: 118mm
White non-woven fabric
Yellow and transparent bag



03-2

外贸系列

(International Range)

二件式造口底盘

Two-piece Skin Barrier

- 产品描述 (Product Description)
- 本产品为回肠、结肠、直肠造口的护理产品。

- 产品接触完好皮肤和肠内腔。

- 非无菌提供。

- This is a medical product for ileum, colon and rectum.

- Product comes into contact with intact skin and intestinal lumen.

- Supplied non-sterile.

■

型号(Model): 7309

剪孔(Cut to fit): 15~57mm

扣件(Flange): 57扣件



■

型号(Model): W7303

剪孔(Cut to fit): 15~75mm

扣件(Flange): 57扣件



■

型号(Model): Y7255

剪孔(Cut to fit): 15~45mm

扣件(Flange): /



■

型号(Model): Y7203

剪孔(Cut to fit): 15~45mm

扣件(Flange): 45扣件



■

型号(Model): Y7401

剪孔(Cut to fit): 15~70mm

扣件(Flange): 70扣件



■ 产品简介 (Product Introduction)

项目	介绍	Item	Introduction
搭配说明	二件式底盘需搭配相同扣件的二件式袋体进行使用。	Matching Instruction	The skin barrier and the same size flange two-piece ostomy bag are used together.
组成结构	二件式底盘由进口水胶体材质、离型纸(膜)和扣件系统组成。	Components	Skin barrier is made of imported hydrocolloid, paper/PET film and flange system.

■ 产品简介 (Product Introduction)

项目	介绍	Item	Introduction
产品特点	独立包装的二件式底盘，粘性好，对皮肤致敏率低，无刺激。扣环与袋体分离的设计，与袋体扣合和拆解时更方便。	Product Features	Individual package of skin barrier, good adhesion, low skin sensitization rate, no irritation. The design that the flange is separated from the bag body is more convenient to fasten and disassemble with the bag body.

03-3



外贸系列

(International Range)

二件式造口袋

Two-piece Ostomy Bag

■ 产品描述 (Product Description)

- 本产品为回肠、结肠、直肠造口的护理产品。

- 产品接触完好皮肤和肠内腔。

- 非无菌提供。

- This is a medical product for ileum, colon and rectum

- Product comes into contact with intact skin and intestinal lumen.

- Supplied non-sterile.

■

型号(Model): Y6201
扣件(Flange): 45扣件





■

型号(Model): Y4308
扣件(Flange): 57扣件





■

型号(Model): W4202
扣件(Flange): 45扣件





■

型号(Model): Y4303B
扣件(Flange): 57扣件





■

型号(Model): Y4402
扣件(Flange): 70扣件





■ 产品简介 (Product Introduction)

项目	介绍	Item	Introduction
搭配说明	二件式袋体需搭配相同扣件的二件式底盘进行使用。	Matching Instruction	A two-piece ostomy bag and the same size flange skin barrier are used together.
产品特点	泌尿袋(Y6201)适用于尿路造口者使用，防逆流设计，有效减少尿液对造口和底盘的腐蚀。	Product Features	Urostomy bag (Y6201) is suitable for urostomy patient, anti-backflow design, effectively reduce the corrosion of urine to stoma and chassis.

■ 产品简介 (Product Introduction)

项目	介绍	Item	Introduction
组成结构	主要由高阻隔膜袋体、无纺布、活性碳片、扣件系统和封口器组成。	Components	The bag is made of EVOH film, non-woven fabric, carbon filter, flange, and closed part.
产品特点	肛肠袋可随时排放内容物并拆解下来清洗再使用，方便实惠。	Product Features	The contents inside the ostomy bags can be discharged at any time and disassembled for cleaning and reuse, which is convenient and affordable.