



CFDA | ISO 9001 | ISO 13485 | CE | FDA



MONITORING ECG ELECTRODES FOR ADULTS

ICU/CCU, SURGERY, EMERGENCY, LONG TERM, Stress Test, Holter, MRI, X-ray

Foam



HS02



HS05




HS28



HS16



HS11



HS06



HS03




HS04



HS30



WS04



WS08



WS28



TS04



TS06



TS02

Cloth



HW04



HW07



HW08




HW09



HR01

Breathable film



HA01



HA07



HA09



HL01



HL02

Prewired Electrode

Monitoring Electrodes offer total versatility at affordable prices. 3 different backings, a standard and an enhanced adhesive the YUECHEN range has an electrode to meet all your particular needs in short or long term monitoring.

Protective Foam Backings of short term monitoring electrodes are impermeable to liquids present in surgical procedures.

Breathable Film or Soft Cloth Backings of long term monitoring electrodes provide patient comfort for long term use.

Solid Gel is not only conductive but also adhesive and can be removed without residue.

Easy to open pouch and a handling aid on all YUECHEN monitoring electrodes facilitate use.


Foam

Code	Gel	Dimensions	Suggested Use
HS02	Solid gel	37*42mm	Monitoring, Short term/Long term
HS05	Solid gel	40*46mm	Monitoring, Short term/Long term
HS28	Solid gel	32*41mm	Monitoring, ICU/CCU, Surgery, Emergency, Long term
HS16	Solid gel	42*45mm	Monitoring, ICU/CCU, Surgery, Emergency, Long term
HS11	Solid gel	50*55mm	Monitoring, ICU/CCU, Surgery, Emergency, Long term
HS06	Solid gel	Φ25.4mm	Monitoring, ICU/CCU, Surgery, Emergency, Long term
HS03	Solid gel	35*44mm	Monitoring, ICU/CCU, Surgery, Emergency, Long term
HS04	Solid gel	Φ50mm	Monitoring, Long term, Stress Test, Holter
HS30	Solid gel	57*43mm	Monitoring, Long term, Stress Test, Holter
WS04	Wet gel	Φ50mm	Monitoring, Long term, Stress Test, Holter
WS08	Wet gel	45*48mm	Monitoring, Long term, Stress Test, Holter
WS28	Wet gel	32*41mm	Monitoring, Long term, Stress Test, Holter
TS04	Solid gel	Φ50mm	X-Ray, MRI, Monitoring during Imaging Procedures
TS06	Solid gel	Φ25.4mm	X-Ray, MRI, Monitoring during Imaging Procedures
TS02	Solid gel	37*42mm	X-Ray, MRI, Monitoring during Imaging Procedures

Cloth, Breathable film, Prewired Electrode

Code	Gel	Dimensions	Suggested Use
HW04	Solid gel	Φ50mm	Monitoring, Short term/Long term
HW07	Solid gel	34*48mm	Monitoring, Short term/Long term
HW08	Solid gel	37*42mm	Monitoring, ICU/CCU, Surgery, Emergency, Long term
HW09	Solid gel	45*48mm	Monitoring, ICU/CCU, Surgery, Emergency, Long term
HR01	Hydrogel	23*34mm	ICU/CCU, Surgery, Emergency, Cath Lab, Radiology
HA01	Solid gel	Φ48mm	Monitoring, ICU/CCU, Surgery, Emergency, Long term, Stress Test, Holter
HA07	Solid gel	34*48mm	Monitoring, ICU/CCU, Surgery, Emergency, Long term, Stress Test, Holter
HA09	Solid gel	45*48mm	Monitoring, ICU/CCU, Surgery, Emergency, Long term
HL01	Solid gel	22*22mm	Monitoring, Long term
HL02	Solid gel	28*35mm	Monitoring, Long term

TAB ELECTRODES / RESTING ECG ELECTRODES



HT01



HT02



HT03



HT04



HT05



HT06

Code	Gel	Dimensions	Suggested Use
HT01	hydrogel	14*30mm	Neonate
HT02	hydrogel	20*30mm	Child/Adult
HT03	hydrogel	26*30mm	Adult
HT04	hydrogel	20*30mm	Child/Adult
HT05	hydrogel	23*34mm	Child/Adult
HT06	hydrogel	28*34mm	Adult

Electrosurgical pencil YUECHEN Brand



HP01 (white, round-button)



HP03 (white, switch-button)



HP02 (blue, round-button)



HP04 (blue, switch-button)



HP05 (blue, oval-button)



HP06 (blue, teflon coated electrode)




HP07 (blue, foot control, Φ6.3mm one pin jack phone plug)






HP08 (blue, foot control, Φ4.0mm one pin plug)



HOLSTER



TIP CLEANER



The YUECHEN Single Use Electrosurgical Pencils are a high quality diathermy pencil. The electrode is fitted with a hexagonal connection, securing the electrode in the required position. It can be used with any electrosurgical generator available.

- Comfortable design and light weight
- Colour marked buttons
- Hexagonal electrode connection
- 3 meter cable, 2.36mm electrode, 70mm blade length (40mm outside of pencil)
- 3 pin US-connector
- Handle length (170mm only for HP05, 153mm the other)

Disposable Neutral Electrodes

The Professional's Choice



HG01



HG02



HG03



HG04



HG05



HG06



HG07



HG11



HG12



HG01-1



HG02-1



HG03-1



HG04-1



HG11-1



HG12-1



ESU Cable




CA01




CA02

Code	Backing	Sort	Dimensions	Size
HG01	Foam	Monopolar	210*100mm	Adult
HG02	Foam	Biopolar	210*100mm	Adult
HG03	Foam	Monopolar	150*90mm	Child
HG04	Foam	Biopolar	150*90mm	Child
HG05	Foam	Monopolar	75*110mm	Neonate
HG06	Foam	Biopolar	75*110mm	Neonate
HG07	Foam	Biopolar	88*128mm	Adult
HG11	Foam	Monopolar	88*128mm	Adult
HG12	Foam	Biopolar	88*128mm	Adult
HG01-1	Foam	Monopolar, REM Plug	210*100mm	Adult
HG02-1	Foam	Biopolar, REM Plug	210*100mm	Adult
HG03-1	Foam	Monopolar, REM Plug	150*90mm	Child
HG04-1	Foam	Biopolar, REM Plug	150*90mm	Child
HG11-1	Foam	Monopolar, REM Plug	88*128mm	Adult
HG12-1	Foam	Biopolar, REM Plug	88*128mm	Adult

TENS Electrodes




TES-B-01




TES-B-02




TES-B-03




TES-B-04



TES-B-06




TES-B-05



TES-B-07


Code	Gel	Dimensions	button
TES-B-01	hydrogel	40*40mm	4.0/5.0
TES-B-02	hydrogel	50*50mm	4.0/5.0
TES-B-03	hydrogel	40*60mm	4.0/5.0
TES-B-04	hydrogel	Φ25mm	4.0/5.0
TES-B-05	hydrogel	Φ32mm	4.0/5.0
TES-B-06	hydrogel	50*90mm	4.0/5.0
TES-B-07	hydrogel	86*50mm	4.0/5.0




TES-L-01

TES-L-02


TES-L-03




TES-L-04



TES-L-06



TES-L-05



TES-L-07

Code	Gel	Dimensions	length of line	Aperture
TES-L-01	hydrogel	40*40mm	100mm or customer requirements	1.0/1.5/2.0/2.5/5.0
TES-L-02	hydrogel	50*50mm	100mm or customer requirements	1.0/1.5/2.0/2.5/5.0
TES-L-03	hydrogel	40*60mm	100mm or customer requirements	1.0/1.5/2.0/2.5/5.0
TES-L-04	hydrogel	Φ25mm	100mm or customer requirements	1.0/1.5/2.0/2.5/5.0
TES-L-05	hydrogel	Φ32mm	100mm or customer requirements	1.0/1.5/2.0/2.5/5.0
TES-L-06	hydrogel	50*90mm	100mm or customer requirements	1.0/1.5/2.0/2.5/5.0
TES-L-07	hydrogel	86*50mm	100mm or customer requirements	1.0/1.5/2.0/2.5/5.0