

(0086571-86822701

Add.: # 701, Building C, Wanda Plaza,No.666,
Hangxing Road,Gongshu District,
Hangzhou, 310015 , Zhejiang Province, China

Web: www.weheremedical.com

Wehere

Hangzhou Wehere Medical Technology Co., Ltd.

Reuscitator



PVC Manual Resuscitator

It is applicable to the occasions of cardiopulmonary resuscitation and respiratory first aid at work. It is especially suitable for asphyxia, dyspnea or the need to increase the oxygen supply. It is easy to use, less pain, less complications, easy to carry, and can be immediately ventilated with or without oxygen source.

Conventional configuration: PVC resuscitation bag, air bag, oxygen tube, PVC anesthesia mask, instructions manual.

Optional accessories: oropharyngeal airway (edge type and central type), opener, PEEP valve.

Packing method: PP box, paper box, PE bag.

- 1. The product uses medical grade PVC. It is non-toxic and non-irritating.
- 2. The mask is made of PVC material, The air cushion contour of the mask fits the human face shape and provides good sealing effect.
- 3.The product specification is complete and is disposable

Item No.	Patient Size
PMR-L	Adult
PMR-M	Pediatric
PMR-S	Infant



Silicone Manual Resuscitator

21

It is especially suitable for asphyxia, dyspnea or the need to increase the oxygen supply. It is easy to use, less pain, less complications, easy to carry, and can be immediately ventilated with or without oxygen source.

Conventional configuration: silicone resuscitation bag, air bag, oxygen tube, silicone anesthesia mask, instructions manual.

Optional 2 with 2 pieces: oropharyngeal airway (edge type and central type), opener, PEEP valve.

Packing method: PP box. Paper box, PE bag.

- 1. The product uses medical grade silicone, which is non-toxic, non-irritating and non-allergic.
- 2. The product specification is complete, and it is reused for 20 times.

Item No.	Patient Size
SMR-L	Adult
SMR-M	Pediatric
SMR-S	Infant

