



(0086571-86822701

Add.: # 701, Building C, Wanda Plaza, No. 666,
Hangxing Road, Gongshu District,
Hangzhou, 310015, Zhejiang Province, China

Web: www.weheremedical.com

Wehere

Hangzhou Wehere Medical Technology Co., Ltd.

Stethoscope





ZY1001 单头听诊器/SINGLE HEAD STETHOSCOPE

铝制单用听头，听头直径46mm，可配铜质耳挂或铝合金耳挂，耳挂配置环保软耳塞。PVC胶管有环保和不环保2种选择，颜色可配（黑、灰、红、绿、蓝等），简易纸盒包装。100pcs/ctn，52x41x28cm，11KGS。

Aluminum single head, chestpiece diameter 46mm, can be equipped with copper ear hanging or aluminum ear hanging, ear hanging environmental protection soft ear tip. PVC tube has two options environmentally and not environmentally, tube color can be available (black, gray, red, green, blue, etc.). Simple box packaging, one piece in one box. 100pcs/ctn 52x41x28cm, 11KGS.



ZY1002 儿童单头听诊器/SINGLE HEAD STETHOSCOPE(PEDIATRIC)

铝制儿童单用听头，听头直径37mm，可配铜质耳挂或铝合金耳挂，耳挂配置环保软耳塞。PVC胶管有环保和不环保2种选择，颜色可配（黑、灰、红、绿、蓝等），简易纸盒包装。100pcs/ctn，52x41x28cm，11KGS。

Child Aluminum single head, chestpiece diameter 37mm, can be equipped with copper ear hanging or aluminum ear hanging, ear hanging environmental protection soft ear tip. PVC tube has two options environmentally and not environmentally, tube color can be available (black, gray, red, green, blue, etc.). Simple box packaging, one piece in one box. 100pcs/ctn 52x41x28cm, 11KGS.



ZY1003 彩色单头听诊器/COLORED SINGLE HEAD STETHOSCOPE

铝制单用听头，听头直径37mm，铝合金耳挂，耳挂配置环保软耳塞。PVC胶管有环保和不环保2种选择，颜色可配（黑、灰、红、绿、蓝等），听头、耳挂、胶管同色，简易纸盒包装。100pcs/ctn，52x41x28cm，11KGS。

Aluminum single head, chestpiece diameter 37mm, can be equipped with copper earSimple box packaging, one piece in one box. 100pcs/ctn 52x41x28cm, 11KGS.hanging or aluminum ear hanging, ear hanging environmental protection soft ear tip. PVC tube has two options environmentally and not environmentally, tube color can be available (black, gray, red, green, blue, etc.). chestpiece, ear hanging, tube are same color. Simple box packaging, one piece in one box. 100pcs/ctn 52x41x28cm, 11KGS.



ZY1004 防冷圈单头听诊器/SINGLE HEAD STETHOSCOPE WITH NON-CHILL RING

铝制单用听头，带防冷圈，听头直径46mm，可配铜质耳挂或铝合金耳挂，耳挂配置环保软耳塞。PVC胶管有环保和不环保2种选择，颜色可配（黑、灰、红、绿、蓝等），简易纸盒包装。100pcs/ctn，52x41x28cm，11KGS。

Aluminum single head with cold proof ring, chestpiece diameter 46mm, can be equipped with copper ear hanging or aluminum ear hanging, ear hanging environmental protection soft ear tip. PVC tube has two options environmentally and not environmentally, tube color can be available (black, gray, red, green, blue, etc.). Simple box packaging, one piece in one box. 100pcs/ctn 52x41x28cm, 11KGS.



ZY1005 碗式听诊器/BOWLES STETHOSCOPE

锌合金特种形状听头，带防冷圈，听头直径46mm，铝合金耳挂，耳挂配置环保软耳塞。PVC胶管颜色可配（黑、灰、红、绿、蓝等），简易纸盒包装。防冷圈胶管同色。100pcs/ctn，52x41x28cm，11KGS。

Zinc alloy special shaped head with cold proof ring, chestpiece diameter is 46mm, can be equipped with copper ear hanging or aluminum ear hanging, ear hanging environmental protection soft ear tip, tube color can be available (black, gray, red, green, blue, etc.). Simple box packaging, cold proof ring color is same as tube. 100pcs/ctn, 52x41x28cm, 11KGS.



ZY1006 ABS电镀单头听诊器/ABS SINGLE HEAD STETHOSCOPE

ABS单头听头，听头直径46毫米，铝质氧化耳挂，配塑料耳塞。PVC连接管，胶管颜色：黑、蓝、红、黄、绿等纸盒包装，100套/箱，52x41x28cm，9公斤

ABS single head, chestpiece diameter 46mm, aluminum ear hanging Plastic ear tip. PVC tube, tube color can be available (black, gray, red, green, blue, etc.). Packing in box, 100sets/box, 52x41x28cm, 9kgs. Simple box packaging, one piece in one box. 100pcs/ctn 52x41x28cm, 11KGS.



ZY-2001 双头听诊器/DAUL HEAD STETHOSCOPE

铝制双用听头，听头直径46mm，可配铜质耳挂或铝合金耳挂，耳挂配置环保软耳塞。PVC胶管有环保和不环保2种选择，颜色可配（黑、灰、红、绿、蓝等）。简易纸盒包装。100pcs/ctn，52x41x28cm，12KGS

Aluminum dual head, chestpiece diameter 46mm, can be equipped with copper ear hanging or aluminum ear hanging, ear hanging environmental protection soft ear tip. PVC tube has two options: environmentally and not environmentally, tube color can be available (black, gray, red, green, blue, etc.) Simple box packaging, one piece in one box. 100pcs/ctn. 52x41x28cm, 11KGS.



ZY-2002 儿童双头听诊器/DAUL HEAD STETHOSCOPE(PEDIATRIC)

铝制双用听头，听头直径37mm，可配铜质耳挂或铝合金耳挂，耳挂配置环保软耳塞。PVC胶管有环保和不环保2种选择，颜色可配（黑、灰、红、绿、蓝等）。简易纸盒包装。100pcs/ctn. 52x41x28cm，12KGS

Aluminum dual head, chestpiece diameter 37mm, can be equipped with copper ear hanging or aluminum ear hanging, ear hanging environmental protection soft ear tip. PVC tube has two options: environmentally and not environmentally, tube color can be available (black, gray, red, green, blue, etc.) Simple box packaging, one piece in one box. 100pcs/ctn. 52x41x28cm, 11KGS.



ZY-2003 彩色双头听诊器/COLORED DAUL HEAD STETHOSCOPE

铝制双用听头，听头直径46mm，铝合金耳挂，耳挂配置环保软耳塞。PVC胶管有环保和不环保2种选择，颜色可配（黑、灰、红、绿、蓝等）。听头、耳挂、胶管同色，简易纸盒包装。100pcs/ctn，52x41x28cm，12KGS

Aluminum dual head, chestpiece diameter 37mm, can be equipped with copper ear hanging or aluminum ear hanging, ear hanging environmental protection soft ear tip. PVC tube has two options: environmentally and not environmentally, tube color can be available (black, gray, red, green, blue, etc.) chestpiece, ear hanging, tube are same color. Simple box packaging, one piece in one box. 100pcs/ctn. 52x41x28cm, 11KGS. 11KGS.



ZY-2004 带防冷圈双头/DUAL HEAD STETHOSCOPE WITH NON-CHILL RING

铝制双用听头，带防冷圈，听头直径46mm，可配铜质耳挂或铝合金耳挂，耳挂配置环保软耳塞。PVC胶管有环保和不环保2种选择，颜色可配（黑、灰、红、绿、蓝等）。简易纸盒包装。100pcs/ctn，52x41x28cm，12KGS

Aluminum dual head with cold proof ring, chestpiece diameter 46mm, can be equipped with copper ear hanging or aluminum ear hanging, ear hanging environmental protection soft ear tip. PVC tube has two options: environmentally and not environmentally, tube color can be available (black, gray, red, green, blue, etc.) Simple box packaging, one piece in one box. 100pcs/ctn, 52x41x28cm, 11KGS.



ZY-3001 多功能听诊器/SPRAGUE RAPPAPORT STETHOSCOPE

全旋转锌合金听头，直径45毫米，铜质镀铬耳挂，22英寸长的双管，配件包含：1对膜片、2对耳塞、三个钟型听头。胶管颜色：黑、红、黄、蓝和透明等，纸盒包装或礼品盒海绵包装。50套/箱，67x31x44cm，19公斤

Full rotation of zinc alloy head, chestpiece diameter 45 mm, copper chrome ear hanging, 22 inch long double tube, accessories include: 1 pairs of diaphragms, 2 pairs of ear tip, three bells. Tube color: black, red, yellow, blue and transparent etc., simple box packaging or gift box sponge packaging. 50 sets / box, 67x31x44cm, 19kgs



ZY-3002 彩色多功能听诊器/SPRAGUE RAPPAPORT STETHOSCOPE

全旋转锌合金听头，直径45毫米，铜质镀铬耳挂，22英寸长的双管，配件包含：1对膜片、2对耳塞、三个钟型听头。听头和胶管同色，胶管颜色：黑、红、黄、蓝和透明等，纸盒包装或礼品盒海绵包装。50套/箱，67x31x44cm，19公斤

Full rotation of zinc alloy head, chestpiece diameter 45 mm, copper chrome ear hanging, 22 inch long double tube, accessories include: 1 pairs of diaphragms, 2 pairs of ear tip, three bells. Chestpiece and tube are same color. Tube color: black, red, yellow, blue and transparent etc., simple box packaging or gift box sponge packaging. 50 sets / box, 67x31x44cm, 19kgs Simple box packaging, one piece in one box. 100pcs/ctn. 52x41x28cm, 11KGS.



ZY-3003 ABS电镀单头听诊器/SPRAGUE RAPPAPORT STETHOSCOPE WITH CLOCK
全旋转锌合金带表听头，直径45毫米，铝质镀铬耳挂，22英寸长的双管，配件包含：1对膜片、2对耳塞、三个钟型听头。听头和胶管同色。胶管颜色：黑、红、黄、蓝和透明管等。纸盒包装或礼品盒海绵包装50套/箱。67x31x44cm, 19公斤

Full rotation of zinc alloy head with a clock/chestpiece diameter 45 mm, copper chrome ear hanging, 22 inch long double tube, accessories include: 1 pair of diaphragms, 2 pairs of ear tip, three bells. Chestpiece and tube are same color. Tube color: black, red, yellow, blue and transparent etc., simple box packaging or gift box sponge packaging. 50 sets / box, 67x31x44cm, 19kgs



ZY-3004 镀金多功能听诊器/GOLD SPRAGUE RAPPAPORT STETHOSCOPE
全旋转镀金锌合金听头，直径45毫米，铜镀金耳挂，22英寸长的双管，配件包含：1对膜片、2对耳塞、三个钟型听头。胶管颜色：黑、红、黄、蓝和透明管等。纸盒包装50套/箱。67x31x44cm, 19公斤

Fully rotation gold plated zinc alloy head with a diameter of 45 mm, copper plated ear hanging, 22 inch double tube, accessories include: 1 pair of diaphragms, 2 pairs of ear tips, and three bell shaped heads. tube colors: black, red, yellow, blue, and transparent tubes, packaged in cardboard boxes 50 sets/box, 67x31x44cm, 19kgs a packaging, one piece in one box. 100pcs/ctn. 52x41x28cm, 11KGS.



ZY-3005 彩虹多功能听诊器/Rainbow SPRAGUE RAPPAPORT stethoscope
全旋转锌合金听头，表面彩虹处理，直径45毫米，黑色耳挂，22英寸长的双管，配件包含：1对膜片、2对耳塞、三个钟型听头。胶管颜色：黑、红、黄、蓝和透明管等。EVA包包装，20套/箱。67x31x44cm, 19公斤

Fully rotation gold plated zinc alloy head with surface rainbow treatment, diameter of 45 mm, black ear hanging, 22 inch double tube, accessories include: 1 pair of diaphragms, 2 pairs of ear tips, and three bell shaped heads. tube colors: black, red, yellow, blue, and transparent tubes, packaged in EVA bag. 20 sets/box, 67x31x44cm, 19kgs



ZY-3006 全黑多功能听诊器/BLACK SPRAGUE RAPPAPORT STETHOSCOPE

铝全旋转锌合金听头，表面黑色处理，直径45毫米，黑色耳挂，22英寸长的双管，配件包含：1对膜片、2对耳塞、三个钟型听头。胶管颜色：黑、红、黄、蓝和透明管等。纸盒包装，50套/箱。67x31x44cm, 19公斤

Fully rotation gold plated zinc alloy head with surface black treatment, diameter of 45 mm, black ear hanging, 22 inch double tube, accessories include: 1 pair of diaphragms, 2 pairs of ear tips, and three bell shaped heads. tube colors: black, red, yellow, blue, and transparent tubes, packaged in EVA bag. 50sets/box, 67x31x44cm, 19kgs



ZY-4001 豪华双头听诊器/DELUXE DUAL HEAD STETHOSCOPE WITH NON-CHILL RING
铝合金带防冷环双用听头，铝制氧化耳挂内簧耳挂，PVC软耳塞。纸盒包装，50套/箱。67x27x44cm, 14公斤

Aluminum alloy with anti-cold ring dual head, copper aluminum plating with yellow ear hanging, PVC soft ear tips, one piece in one box 50 sets / box, 67x27x44cm, 14kgs



ZY-4002 全彩豪华双头听诊器/COLORED DELUXE DUAL HEAD STETHOSCOPE
铝合金带防冷环双用听头，铝制氧化耳挂内簧耳挂，PVC软耳塞。听头、耳挂、胶管同色。纸盒包装，50套/箱。67x27x44cm, 14公斤

Aluminum alloy with anti-cold ring dual head, copper aluminum plating with yellow ear hanging, PVC soft ear tips. Whole stethoscope was same color, one piece in one box 50 sets / box, 67x27x44cm, 14kgs



ZY-4003B 彩虹豪华单用听诊器/Rainbow luxury single head stethoscope

锌合金听头，表面彩虹处理，豪华高档感觉，表面黑色处理内簧耳挂，PVC软耳塞，海绵礼品纸盒包装。25套/箱，84x34x25cm，14公斤

Zinc alloy head, Surface rainbow treatment, Luxury and atmospheric, with black treated inner spring ear hanging and PVC soft ear tips, Sponge gift box packaging. 25 sets / box, 84x34x25cm, 14kgs



ZY-4003 豪华单头听诊器/DELUXE SINGLE HEAD STETHOSCOPE

锌合金听头，表面沙钱处理，豪华高档感觉，铝合金亚光氧化处理内簧耳挂。

Zinc alloy head, the surface is sand nickel treated, luxury surface, the ear hanging are copper with sand nickel surface treated



ZY-4003A 豪华全黑单用听诊器/DELUXE BLACK SINGLE HEAD STETHOSCOPE

锌合金听头，表面亚光黑处理，豪华高档感觉，表面黑处理内簧耳挂，PVC软耳塞，海绵礼品纸盒包装。25套/箱，84x34x25cm，14公斤

Zinc alloy head, Surface matte black treatment, Luxury and atmospheric, with black treated inner spring ear hanging and PVC soft ear tips, Sponge gift box packaging.



ZY-4004 豪华黑色听诊器/DELUXE BLACK STETHOSCOPE

锌合金听头，新型设计，表面精致黑色处理，配高档变频膜片，内置弹簧耳挂，PVC软耳塞或硅胶软耳塞，海绵礼品纸盒包装。25套/箱，84x34x25cm，14公斤

Zinc alloy head, new design, the surface is black treatment, with high-grade frequency diaphragm, built-in spring ear hanging, PVC soft ear tips or soft silicone ear tips, sponge gift box packaging. 25 sets / box, 84x34x25cm, 14kgs



ZY-4005 一次性听诊器/DISPOSABLE STETHOSCOPE

听头、耳挂、胶管全部塑料听诊器，尼龙袋包装，200个/箱 54x42x27cm，11公斤

Chestpiece, ear hanging, tube are all plastic stethoscope. Nylon bag packaging. 200 / box, 54x42x27cm, 11 kgs



ZY-4006 木质胎音听诊器/WOODEN FETAL STETHOSCOPE

桦木制造，专业产科使用纸盒包装，100套/箱，52x27x35cm，8公斤

Made of Birch, professional obstetric using one piece in one box. 100 sets / box, 52x27x35cm, 8 kgs



ZY-4007 塑料胎音听诊器/Plastic FETAL STETHOSCOPE

塑料制造，专业产科使用纸盒包装，100套/箱，52x27x35cm，8公斤

Made of Plastic, professional obstetric use one piece in one box. 100 sets / box, 52x27x35cm, 8 kgs



ZY-4008 亚克力听诊器/Acrylic stethoscope

亚克力听头。铜质镀铬或不锈钢耳挂。PVC软耳塞。海绵礼品包装。25套/箱。84x34x25cm。13公斤

Acrylic head, copper chrome or stainless steel ear hanging. PVC soft ear tips, sponge gift box packaging 25 sets / box. 84x34x25cm, 13kgs



ZY-4009 教学双头听诊器/TEACHING DUAL HEAD STETHOSCOPE

双耳挂。铝制双头听诊器。带加长型PVC连接管。供教学使用。

Double ear hangings, aluminum double head stethoscope with extended PVC connection tube for teaching use.



ZY-5001 心脏科不锈钢听诊器/Cardiology Stainless Steel Stethoscope

深锥型不锈钢听头。适合心脏听诊。配内簧不锈钢耳挂。特种壁厚双腔管。配柔软硅胶耳塞。礼品海绵和包装。25套/箱。84x34x25cm。16公斤

Deep-cone stainless steel head. for heart auscultation, with internal spring stainless steel ear hanging, special wall thickness double tube, with soft silicone ear tips. Gift sponge box packaging. 25 sets / box. 84x34x25cm, 16kgs



ZY-5002 双面膜心脏科听诊器/DOUBLE DIAPHRAGM Cardiology Stainless Steel Stethoscope

双膜片不锈钢听头。适合心脏听诊。配内簧不锈钢耳挂。特种壁厚双腔管。配柔软硅胶耳塞。礼品海绵和包装。25套/箱。84x34x25cm。16公斤

Double diaphragm stainless steel head, suitable for heart auscultation, with internal spring stainless steel ear hanging, special wall thickness double tube, with soft silicone ear tips. Gift sponge box packaging. 25 sets / box. 84x34x25cm, 16kgs



ZY-5003 豪华心脏科听诊器/DELUXE Cardiology Stainless Steel Stethoscope

双膜片不锈钢听头。膜片为新型设计连接膜片。密封性极好。配内簧不锈钢耳挂。特种壁厚双腔管。配柔软硅胶耳塞。礼品海绵和包装。25套/箱。84x34x25cm。16公斤

Double diaphragm stainless steel head. new design Siamese diaphragm, excellent sealing, with internal spring stainless steel ear hanging, special wall thickness double tube, with soft silicone ear tips. Gift sponge box packaging same as tube. 100pcs/ctn. 52x41x28cm, 11KGS 25 sets/ box. 84x34x25cm, 16kgs



ZY-5004 成人不锈钢听诊器/Stainless Steel Stethoscope

经典不锈钢双用听头。内簧不锈钢耳挂。PVC胶管。柔软硅胶耳塞。礼品海绵和包装。25套/箱。84x34x25cm。13公斤

Classic stainless steel dual head, inner spring stainless steel ear hanging. PVC tube, soft silicone ear tips. Gift sponge box packaging. 25 sets / box. 84x34x25cm, 13kgs Simple box packaging one piece in one box. 100pcs/ctn. 52x41x28cm, 11KGS.



ZY-5005 豪华不锈钢听诊器/DELUXE Stainless Steel Stethoscope

新型设计双面膜片不锈钢听头，膜片为连体膜片，内簧不锈钢耳挂，PVC胶管，柔软硅胶耳塞。礼品海绵和包装。25套/箱，84x34x25cm，13公斤

New design double diaphragm stainless steel head, new design Siamese diaphragm, the spring stainless steel ear hanging, PVC tube, soft silicone ear tips. Gift sponge box packaging. 25 sets/ box, 84x34x25cm, 13kgs



ZY-5005A 全黑豪华成人不锈钢听诊器/BLACK Cross-shaped membrane lock Stainless Steel Stethoscope

新型设计双面膜片不锈钢听头，听头表面亚光黑色处理，膜片为新设计连体膜片，密封性极好合金黑色耳挂，PVC胶管。礼品海绵盒包装。25套/箱，84x34x25cm，13公斤

New design double-sided diaphragm stainless steel head. The head surface is matte black treatment. The diaphragm is a newly designed conjoined diaphragm with sealing aluminum alloy black ear hanging and PVC tube. Gift Box Packaging. 25 sets/box, 84x34x25cm, 13kgs 25 sets / box, 84x34x25cm, 16kgs



ZY-5006 变频膜片不锈钢听诊器/Cross-shaped membrane lock Stainless Steel Stethoscope

经典不锈钢双用听头，听头配高档变频膜片，内簧不锈钢耳挂，PVC胶管，柔软硅胶耳塞。礼品海绵和包装。25套/箱，84x34x25cm，13公斤

Classic stainless steel dual head, chestpiece with high-end frequency diaphragm, spring stainless steel ear hanging, PVC tube, soft silicone ear tips. Gift sponge box packaging 25 sets / box, 84x34x25cm, 13kgs



ZY-5007 儿童不锈钢听诊器/Stainless Steel Stethoscope (PEDIATRIC)

经典儿童不锈钢双用听头，内簧不锈钢耳挂，PVC胶管，柔软硅胶耳塞。礼品海绵和包装。25套/箱，84x34x25cm，12公斤

Classic stainless steel dual head (Pediatric), chestpiece with high-end frequency diaphragm, spring stainless steel ear hanging, PVC tube, soft silicone ear tips. Gift sponge box packaging 25 sets/ box, 84x34x25cm, 12kgs



ZY-5008 特种铝合金单头听诊器/Special aluminum alloy single head stethoscope

铝合金特种设计听头，防冷圈，变频膜片，铝合金内簧耳挂，硅胶软耳塞，礼品海绵和包装。25套/箱，84x34x25cm，13公斤

Special design aluminum alloy head, cold proof ring, frequency conversion diaphragm, aluminum alloy inner spring ear hanging, silicone soft ear tips, gift box packaging. 25 sets/box, 84x34x25cm, 13kgs 25 sets / box, 84x34x25cm, 13kgs



ZY-5009 全黑成人不锈钢听诊器/BLACK Stainless Steel Stethoscope

经典不锈钢双用听头，听头表面亚光黑色处理，内簧黑色耳挂。礼品海绵盒包装。25套/箱，84x34x25cm，13公斤

Classic stainless steel dual head stethoscope, with a matte black surface treatment and a black inner spring ear hanging. Gift box packaging. 25 sets/box, 84x34x25cm, 13kgs use.



ZY-5010 A型单头听诊器/A-type single head stethoscope

锌合金电镀头。铝合金电镀内簧耳挂。硅胶软耳塞。A型壁厚双腔管。变频膜片。礼品海绵盒包装。25套/箱。84x34x25cm, 13公斤

Zinc alloy electroplating head, aluminum alloy electroplating inner spring ear hanging, silicone soft ear tips, A-type double tube, frequency conversion diaphragm. Gift box packaging. 25 sets/box, 84x34x25cm, 13kgs



ZY-5011 彩虹头双面不锈钢听诊器/Rainbow head double-sided diaphragm stainless steel stethoscope

新型设计双面膜片不锈钢听头。听头彩虹处理。膜片为新设计连体膜片。密封性极好。耳挂玫瑰金表面处理。PVC胶管。礼品海绵盒包装。25套/箱。84x34x25cm, 13公斤

New design double-sided diaphragm stainless steel head with rainbow treatment, new design conjoined diaphragm, excellent sealing, rose gold surface treatment for ear hanging, PVC tube. Gift Box Packaging 25 sets/box, 84x34x25cm, 13kgs



ZY-5012 咖啡色听诊器/Full Coffee stethoscope

新型设计双面膜片不锈钢听头。听头咖啡色处理。膜片为新设计连体膜片。密封性极好。耳挂咖啡色表面处理。PVC胶管。礼品海绵盒包装 25套/箱。84x34x25cm, 13公斤

New design double-sided diaphragm stainless steel head with coffee color treatment, new design conjoined diaphragm, excellent sealing, coffee color surface treatment for ear hanging, PVC tube. Gift Box Packaging 25 sets/box, 84x34x25cm, 13kgs



ZY-5013 彩虹头黑色耳挂听诊器/Rainbow head black stethoscope

新型设计双面膜片不锈钢听头。听头彩虹处理。膜片为新设计连体膜片。密封性极好。耳挂黑色表面处理。黑色PVC胶管。礼品海绵盒包装 25套/箱。84x34x25cm, 13公斤

New design double-sided diaphragm stainless steel head, with rainbow treatment, new design conjoined diaphragm, excellent sealing, black surface treatment for ear hanging, black PVC tube. Gift Box Packaging 25 sets/box, 84x34x25cm, 13kgs



ZY5014 全铜双头听诊器/Copper Made Dual head stethoscope

铜电镀带防冷环双用听头。铜质电镀内簧耳挂。PVC软耳塞。纸盒包装 50套/箱。67x27x44cm, 14公斤

Copper plated dual head with cold proof ring, copper plated inner spring ear hanging, PVC soft ear tips, box packaging 50 sets/box, 67x27x44cm, 14kgs



ZY5014A 铜彩虹双头听诊器/Copper Made rainbow Dual head stethoscope

铜电镀带防冷环双用听头。听头表面彩虹处理。铜质电镀内簧耳挂。PVC软耳塞。纸盒包装 25套/箱。84x34x25cm, 14公斤

Copper plated dual head with cold proof ring, rainbow treated surface head, copper plated inner spring ear hanging, PVC soft ear tips, box packaging 25 sets/box, 84x34x25cm, 14kgs



ZY-6001 叩诊锤/NEUROLOGICAL HAMMER

锌合金圆头。塑料柄

Zinc alloy round head, plastic handle



ZY-6002 全黑叩诊锤/Black NEUROLOGICAL HAMMER

黑色锌合金圆头。塑料柄

Black zinc alloy round head, plastic handle




ZY-6003 叩诊锤/NEUROLOGICAL HAMMER

锌合金圆头，塑料柄 Zinc alloy round head, plastic handle.



ZY-6004 叩诊锤/NEUROLOGICAL HAMMER

锌合金圆头，塑料柄 Zinc alloy round head, plastic handle.



ZY-6005 六角柄叩诊锤/Hexagonal NEUROLOGICAL HAMMER

锌合金圆头电镀头，电镀金属六角柄

Zinc alloy round electroplating head, electroplating metal hexagonal handle



ZY-6006 黑色六角柄叩诊锤/Black Hexagonal NEUROLOGICAL HAMMER

黑色锌合金圆头，黑色金属六角柄

Black zinc alloy round head, black metal hexagonal handle



ZY-8001 彩管彩虹双面不锈钢听诊器/Colorful tube rainbow double-sided diaphragm stainless steel stethoscope

新型设计双面膜片不锈钢听头，听头彩虹处理，膜片为新设计连体膜片，密封性极好，耳挂彩虹表面处理，新型花色胶管。礼品海绵盒包装。25套/箱，84x34x25cm，13公斤

New design double-sided diaphragm stainless steel head with rainbow treatment, new design conjoined diaphragm, excellent sealing, rainbow surface treatment for ear hanging, new colorful tube. Gift Box Packaging. 25 sets/box, 84x34x25cm, 13kgs



ZY-8002 彩管咖啡色双面不锈钢听诊器/Colorful tube coffee double-sided diaphragm stainless steel stethoscope

新型设计双面膜片不锈钢听头，听头咖啡色处理，膜片为新设计连体膜片，密封性极好，耳挂咖啡色表面处理，新型花色胶管。礼品海绵盒包装。25套/箱，84x34x25cm，13公斤

New design double-sided diaphragm stainless steel head with coffee color treatment, new design conjoined diaphragm, excellent sealing, coffee color surface treatment for ear hanging, new colorful tube. Gift Box Packaging. 25 sets/box, 84x34x25cm, 13kgs



ZY-8003 彩管电镀豪华单听/Colorful tube electroplating luxury single head stethoscope

锌合金电镀头，铝合金电镀内簧耳挂，硅胶软耳套，新型彩色胶管，变频膜片。礼品海绵盒包装。25套/箱，84x34x25cm，13公斤

Zinc alloy electroplating head, aluminum alloy electroplating inner spring ear hanging, silicone soft ear tips, new colorful tube, frequency conversion diaphragm. Gift box packaging. 25 sets/box, 84x34x25cm, 13kgs



ZY-8004 彩管全黑成人不锈钢听诊器/Colorful tube black stainless steel stethoscope

经典不锈钢双用听头，听头表面亚光黑色处理，内簧黑色不锈钢耳挂。礼品海绵盒包装。25套/箱，84x34x25cm，13公斤

Classic stainless steel dual head with matte black surface treatment and black stainless steel inner spring ear hanging. Gift box packaging. 25 sets/box, 84x34x25cm, 13kgs